
A Computer Program for Processing Impedance Cardiographic Data: Improving Accuracy Through User- Interactive Software

Patricia S. Cowings, Karen Naifeh,
and Chet Thrasher

(NASA-TM-101020) A COMPUTER PROGRAM FOR
PROCESSING IMPEDANCE CARDIOGRAPHIC DATA:
IMPROVING ACCURACY THROUGH USER-INTERACTIVE
SOFTWARE (NASA) DC 7 CSCI 06P

N89-12192

Unclas
G3/52 0174649

October 1988



National Aeronautics and
Space Administration



A Computer Program for Processing Impedance Cardiographic Data: Improving Accuracy Through User- Interactive Software

Patricia S. Cowings, Ames Research Center, Moffett Field, California

Karen Naifeh

Chet Thrasher, University of California, San Francisco, California

October 1988



National Aeronautics and
Space Administration

Ames Research Center
Moffett Field, California 94035

A COMPUTER PROGRAM FOR PROCESSING IMPEDANCE CARDIOGRAPHIC DATA:

IMPROVING ACCURACY THROUGH USER-INTERACTIVE SOFTWARE

Patricia S. Cowings, Karen Naifeh,* and Chet Thrasher*

Ames Research Center

SUMMARY

This report contains the source code and documentation for a computer program used to process impedance cardiography data. The cardiodynamic measures derived from impedance cardiography are ventricular stroke volume, cardiac output, cardiac index and Heather index. The program digitizes data collected from the Minnesota Impedance Cardiograph, electrocardiography (ECG), and respiratory cycles and then stores these data on hard disk. It computes the cardiodynamic functions using interactive graphics and stores the means and standard deviations of each 15-second data epoch on floppy disk. This software was designed on a Digital PRO380 microcomputer and used version 2.0 of P/OS, with (minimally) a 4-channel 16 bit analog/digital (A/D) converter. Applications software is written in Fortran 77, and uses Digital's Pro-tool Kit Real Time Interface Library (PRTIL), CORE Graphic Library (CGL), and laboratory routines. Source code can be readily modified to accommodate alternative detection, A/D conversion and interactive graphics. The object code utilizing overlays and multitasking has a maximum of 50 Kbytes.

INTRODUCTION

The Psychophysiological Research Laboratory of the Neurosciences Branch at NASA Ames Research Center has been engaged in a series of ground-based investigations of human autonomic responses to motion sickness. With the final goal of developing a treatment for the motion sickness-like symptoms which affect astronauts exposed to the microgravity environment of space, our group uses noninvasive electrophysiological measures to document changes in physiological activity levels in different subject populations. In the course of this research, we have found that measures of cardiovascular function are very sensitive indices of the malaise levels experienced by test participants. In previous studies (ref. 1), it was observed that high-susceptibility subjects tended to produce more labile and larger magnitude changes in heart rate and blood volume pulse (a relative measure of peripheral resistance), when exposed to motion sickness stressors than low-susceptibility subjects. We decided to investigate this result further.

*University of California, San Francisco, CA.

Although electrocardiography (ECG) provides excellent information regarding the electrical function of the heart muscle, it gives no definitive information regarding physical function. Heart rate alone is insufficient for determining the level of sympathetic influence on cardiac muscle. An increase in heart rate can be caused by either an increase in sympathetic tone or a decrease in parasympathic activation. Finger pulse volume, although easy to monitor, is not as reliable an index of sympathetic tone of the vasomusculature as total peripheral resistance, which is derived from both blood pressure and cardiac output (ref. 2).

We needed a noninvasive, atraumatic means of examining human cardiodynamics in a motion-sickness inducing environment. Our assessment of impedance cardiography is that it provides relatively accurate data on a wide range of dynamic cardiovascular variables and appears even to be sensitive to rapid changes in these variables. Impedance cardiography techniques enable the measurement of myocardial contractility, (heather index) which is directly related to the level of sympathetic innervation of the heart. Even when there is no apparent change in heart rate, changes in stroke volume or myocardial contractility may be occurring which reflect significant differences in autonomic function (refs. 2,3). Other valuable information obtained from this device include cardiac output and central volume (transthoracic impedance).

Miller and Horvath (ref. 3) describe the advantages and drawbacks of impedance cardiography and compare the accuracy of this technique to other invasive and noninvasive measures of cardiodynamics. The principal disadvantage is that erroneous results may be produced if the user is unfamiliar with the effects of artifact on the computation of specific data epochs and is not careful to select cardiac cycles (dZ/dT peaks) that occur only during the expiratory pause between breaths. The Minnesota Impedance Cardiograph, manufactured by Instrumentation for Medicine, Inc., Greenwich, Conn., is available with a (firmware only) microprocessor which is designed to automate this process. However, we found that with the microprocessor, one can obtain data from one heartbeat at a minimum of 8-sec intervals. Thus much data are lost from the intervening beats, especially when rapid changes are occurring in the cardiovascular system.

This paper contains the source code and operator's instructions for a user-interactive program written in support of our research. The program provides more accurate calculations of cardiac parameters based on a greater quantity of data. The sampling rate for digitizing data is 200 samples/sec; the data can be sampled at twice real-time speed from analog tape. The collected data consist of the first derivative of pulsatile thoracic impedance change (dZ/dt), basal impedance, electrocardiogram, and respiration waveforms. During an interactive graphics session, data for dZ/dt , ECG and chest respiration are displayed on a monitor in 15-sec epochs. The user can select from the screen the dZ/dt peaks that are used to calculate cardiodynamic functions. The program computes and stores the means and standard deviations for the cardiac measures for each 15-sec epoch.

CALCULATIONS PERFORMED:

$$\text{Stroke Volume} = \frac{PL^2 T(dZ/dt)_{\min}}{Z_0^2}$$

$$P = 53.2 e^{(0.002)(\text{Hematocrit})}$$

Expressed in ohms-cm

$e = 2.7183$ the natural exponent

Hematocrit expressed in units of percent

L = distance between two inner electrodes in cm

Z_0 = Transthoracic Impedance expressed in ohms

$$\text{Cardiac Output} = (\text{Stroke Volume})(\text{Pulse rate})/1000$$

Expressed in liters per min

Pulse rate is expressed in beats per min

$$\text{Cardiac Index} = \text{Cardiac Output}/\text{Body Surface Area}$$

Expressed in liters per min

$$\text{Body surface area} = H \times W \times 0.007184$$

Expressed in m^2

H = height in cm

W = weight in Kg

Heather Index = $(dZ/dt)_{\min}/R-Z$

Expressed in ohms per sec²

R-Z = interval between the R-wave of the ECG and
the peak of dZ/dt expressed in sec

REQUIRED HARDWARE AND SUBJECT INFORMATION

Hardware for this research includes a Minnesota Impedance Cardiograph, a respiration transducer (e.g., a piezoelectric or mercury strain gage) and a preamplifier capable of producing an analog output signal of respiratory responses. electrocardiography equipment (either the Minnesota Impedance cardiograph for direct measurement or an ECG amplifier for external measurement), a Digital PRO380 microcomputer which uses version 2.0 of P/OS, and a four-channel A/D converter (16-bit). Applications software, written in Fortran 77, uses Digital's Pro-tool Kit Real Time Interface Library (PRTIL), CORE Graphic Library (CGL), and laboratory routines. Source code can be readily modified to accommodate alternative routines for peak detection, A/D conversion and interactive graphics. The object code uses overlays and multitasking and has a maximum of 50 Kbytes.

Additional subject information required to implement this software is (a) hematocrit count; (b) weight in Kg; (c) height in cm; and (d) the distance between the two inner impedance cardiography electrodes (tapes), measured both in front and back and then averaged.

This software is installed on a DEC PRO 380 computer by following applications installation instructions in Professional Developer's Toolkit Reference Manual, Chapters 3 and 6.

ORIGINAL PAGE IS
OF POOR QUALITY

OPERATOR INSTRUCTIONS

Three options will appear on the screen:

```
*****
*
*  SELECT FROM THE FOLLOWING
*
*  1. DIGITIZE A DATA FILE FROM TAPE
*
*  2. DATA REDUCTION OF DIGITIZED FILE
*
*  3. EXIT
*
*****
```

[ENTER 1, 2, OR 3]:

Procedure for digitizing a tape data file--SELECTION 1. If the user selects option 1, the following displays appear on the screen:

```
*****
                        DIGITIZING A DATA FILE
*****
```

PLEASE ENTER RUN ID FOR THIS FILE:

1. User enters an ID of 6 digits.

PLEASE ENTER RUN TIME (IN SECS.) FOR TAPE FILE:

2. Determine maximum duration of data file in seconds (e.g. 30 minutes=1800 seconds); ADD 60 seconds to run length (e.g. , 1860). ENTER THIS NUMBER. NOTE: an additional 60 seconds is added to the file length to accommodate acquisition of calibration data.

HOW FAST TO SAMPLE THE DATA ON TAPE?

ENTER 1 FOR REAL-TIME

ENTER 2 FOR TWICE REAL TIME

[ENTER 1 OR 2]:

3. The option of digitizing data at twice real-time is determined by the memory capacity (512 K RAM) and the I/O response time (less than 0.025 sec) of the system implementing this program. This I/O response time is required to store four-channels of data, in two-byte segments. After entering 1 or 2, the screen displays:

HIT "S" TO START, "P" TO PAUSE, OR "A" TO ABORT:

FOLLOWED BY A [RETURN]

4. Position analog tape to beginning of run, then press the keys "S" and "Return" to begin computer acquisition of high calibration data. Data acquisition may be paused at any time by pressing the key "P" and will not continue until the user presses the "RESUME" key. Acquisition may be aborted at any time by pressing the key "A", at which time the program prompts to the user to either SAVE or DELETE the created digitized file. If the user has chosen to start data acquisition (i.e., have pressed the key "S"), the screen will display a flashing message and three columns will scroll continuously until acquisition ends. The first column indicates the active buffer (switches between buffers 1 and 2). The second column indicates error status (e.g., if A/D is turned off). No error is indicated by "0". The third column displays the data second currently being acquired.

ORIGINAL TESTED
OF POOR QUALITY

* DATA ACQUISITION IN PROGRESS *

1	0	1
1	0	1
1	0	1
1	0	1
1	0	1
2	0	2
2	0	2
2	0	2
:	:	:
1	0	30
1	0	30
1	0	30

IMPACQ -- PAUSE - ADVANCE TAPE TO START OF LOW CALS NOW,

THEN PRESS RESUME TO CONTINUE

5. After 30 seconds of high calibration data, the program will automatically pause. Advance tape to start of low calibration data. Press the key "RESUME". Again, the display indicates that data acquisition is in progress.

6. After 30 seconds of low calibration data, the program will automatically pause. When tape has been positioned at start of the data session, computer acquisition is started by pressing the key "RESUME".

7. When the total number of seconds defining this file have elapsed the program screen display will direct the user to: SAVE the file; or DELETE the created digitized file. The screen displays:

SELECT ONE OF THE FOLLOWING [ENTER 1 OR 2]

1 TO SAVE DIGITIZED FILE

2 TO DELETE DIGITIZED FILE

8. If the user has chosen to SAVE the file, it is written to hard disk, and the user may now proceed with data reduction to calculate stroke volume, cardiac output, cardiac index and Heather index. The original menu is displayed:

```
*****
*
*  SELECT FROM THE FOLLOWING
*
*  1. DIGITIZE A DATA FILE FROM TAPE
*
*  2. DATA REDUCTION OF DIGITIZED FILE
*
*  3. EXIT
*
*****
```

[ENTER 1, 2, OR 3]:

Procedure for data reduction--SELECTION 2

1. The screen will request input from the user as follows:

PHASE 1 -- KEY FIELD ENTRY

ENTER RUN ID FOR THIS FILE: (6 digit file name used when digitizing).

ENTER LENGTH (CM) BETWEEN INNER ELECTRODES: (e.g., 26.75)

ENTER HEMATOCRIT COUNT: (e.g., 40)

ENTER HEIGHT (CM): (e.g., 183)

ENTER WEIGHT (KG): (e.g., 94.5)

PHASE 2 -- CALIBRATION VALUES

2. The display will show the calibration levels set internal to the program which the user may choose to modify. These default values are:

	High -----	Low -----
BASE IMPEDANCE (OHMS)	25.0000	0.0000
dZ/dT (OHMS)	1.0000	0.0000

WOULD YOU LIKE TO CHANGE THESE VALUES (Y/N)

NOTE: If the user enters "N", the screen will display:

PHASE 3 -- CALIBRATION ACQUISITION

CALIBRATION OF BASELINE IMPEDANCE

CHANNEL -----	SIGNAL -----	A/DHIGH -----
1	IMPEDANCE	5400

DO YOU WANT TO RE-RUN HIGH CALS [Y/N]?

NOTE: If the user enters "Y" then the the program will return to the original menu and user must redigitize calibration data from tape. The condition that would require a "Y" response here is obtaining an A/DHIGH value which does not correspond to the voltage out of the Minnesota Cardiograph's internal "Hi CAL" setting. In this example, an A/D value of 5400 equals 0.8 volts. For additional information on determining the ratio of A/D values to voltage, refer to the: Pro/Tool Kit Real-Time Interface Library (PRTIL) User's Guide, Chapter 7, p. 12.

If the user enters "N" then the screen displays

CHANNEL -----	SIGNAL -----	A/DLOW -----
1	IMPEDANCE	433

DO YOU WANT TO RE-RUN LOW CALS [Y/N]?

NOTE: If the user enters "Y" then the the program will return to the original menu and user must redigitize calibration data from tape. If the user enters "N" then the screen displays:

CHANNEL	SIGNAL	ADLOW	ADHI	LOWCAL	HICAL
-----	-----	-----	-----	-----	-----
1	IMPEDANCE	433	5400	0.0000	25.0000

WOULD YOU LIKE TO REPEAT THESE CALIBRATIONS [Y/N]?

NOTE: If the user enters "Y" then the the program will return to the original menu and user must redigitize calibration data from tape. If the user enters "N" then the screen displays

CALIBRATION OF dz/dt SIGNAL

CHANNEL	SIGNAL	ADLOW	ADHI	LOWCAL	HICAL
-----	-----	-----	-----	-----	-----
2	dz/dt	491	1495	0.0000	1.0000

SLOPE 9.960159E-04 INTERCEPT -0.4890438

WOULD YOU LIKE TO REPEAT CALIBRATIONS [Y/N]?

NOTE: If the user enters "Y" then the program will return to the original menu and user must redigitize calibration data from tape. If the user enters "N" and this is the first time these data are to be analyzed, the program will display an informative error message: -FILE DOES NOT EXIST-. It will then create an interim file and proceed directly to DATA REDUCTION AND EDITING. If this is an editing session, the screen will display a list of all 15-second data epochs from which the user has selected at least one waveform for processing. If no waveform was chosen by the user for a particular epoch, that epoch is "ZEROED" (i.e., the epoch number is not listed below and will not be used). A 10-min file where all epochs contain selected data is displayed as follows:

FROM THE RECORD, THE FOLLOWING EPOCHS ARE NONZERO:

1	2	3	4	5	6	7	8	10	11
12	13	14	15	16	17	18	19	20	21
22	23	24	25	26	27	28	29	30	31
32	33	34	35	36	37	38	39	40	

PRESS [RETURN] TO CONTINUE

ORIGINAL PAGE IS
OF POOR QUALITY

When the [RETURN] key is pressed, the screen will display:

DATA REDUCTION AND EDITING

ENTER STARTING TIME CODE FOR DATA ON TAPE (e.g., enter 300 for 3:00).

ENTER EPOCH FOR START OF DATA REDUCTION: (e.g., 1)

ENTER LAST DATA EPOCH: (e.g., 40)

The time code entry is recordkeeping information which is not used by the program. It is here for the convenience of the user so that he may keep track of "where" in a particular tape file data reduction was started. If this is an editing session, the user may choose to enter only those epochs needing correction, i.e., new waveforms will be selected. At this point, the screen will display the data of the first 15-sec epoch selected:




Figure 1.- Example of a screen display showing a 15-sec data epoch.

The first waveform displayed in this figure is the dZ/dt . The horizontal lines through this waveform represent high cal level (upper line) and the low cal level (lower line). The two horizontal dots beneath each dZ/dt waveform mark the point where the rising (major peak) dZ/dt waveform first crosses the low-cal line and the maximum trough (downward spike) following this peak. The distance between these two dots (for a single waveform) is the ventricular ejection time.

The numerics displayed beneath the dZ/dt identify this epoch number and time code at the start of this epoch. In this figure, the display 1 300 means epoch 1, time code 3:00.

The second waveform is the electrocardiogram (ECG) and the horizontal dots beneath the ECG mark the peak of the R-wave. The third waveform shows respiratory cycles (Note that this waveform is inverted, such that the top of each cycle displayed actually represents the expiratory pause between breaths). Vertical dashed lines denote the dZ/dt peaks and their relation to ECG and respiration signals. Each dZ/dt peak (vertical dashed lines) is numbered at the bottom of the screen.

While this display is on the screen, the user is prompted with questions (displayed in upper left corner). The first question is:

SKIP THIS EPOCH? [ENTER Y/N]:

If the user does not want to have this data epoch represented in the stored output file, entering "Y" would cause the program to display the second 15-sec data epoch. No data from the first epoch are stored, and the second epoch is renumbered as epoch 1. If the user enters "N" in response to this question, the next query displayed is

REDO ANALYSIS? [ENTER Y/N]:

One condition that would make the user choose to redo the analysis is that the dZ/dt waveform is significantly distorted by artifact (e.g., electrical noise or movement) and therefore the horizontal dots beneath the waveform are incorrectly positioned (i.e., ventricular ejection time is incorrectly labeled). Answering "Y" to this question results in the prompt

ENTER GATE AND DIFFERENCE FACTORS SEPARATED WITH A COMMA:

NOTE: The GATE and DIFFERENCE factors are indigenous to the peak detection subroutine used in this program. They are used for determining the peak and delta-t (ventricular ejection time) values of the dZ/dt waveform. Both factors provide a criterion for distinguishing between high frequency noise and significant trends in the waveform. The GATE factor represents the minimum number of consecutive data points (A/D values) needed to establish a trend. The DIFFERENCE factor represents the minimum difference in magnitude between data points that is needed to establish a real increase in the trend. If the current maximum A/D value is greater than the current maximum plus the minimum DIFFERENCE (current minus previous maximum A/D

value), then the counter is incremented. When the counter is equal to the GATE, a trend toward significant peak is established. The default value for GATE is 2 and for DIFFERENCE is 100. If the user feels that the points marking delta-t are "too close together", then a higher DIFFERENCE factor is entered (e.g., 2,200 or 2,300). If the delta-t marks are "too far apart", then a lower DIFFERENCE factor is entered (e.g., 2,50 or 2,75). The GATE factor is rarely modified. After new GATE or DIFFERENCE factors have been entered, the screen blanks, and this epoch is again displayed with the dots marking delta-t repositioned. This process can be repeated until the user is satisfied that most dZ/dt waveforms in this epoch have delta-t correctly marked. If the user enters "N", the next query displayed is

OK TO CALCULATE HEATHER INDEX? [ENTER Y/N]:

If the user enters "N", the output file would contain zeroes for the heather index. The condition that would make the user choose to zero the heather index is when the ECG waveform is so distorted by artifact that the peak of the R-wave cannot be detected (i.e., marked by a dot) by the program. If the user enters "Y", the output file will contain the calculated heather index for those waveform peaks selected below. The next instruction from the program is

LIST SELECTED PEAKS, SEPARATE WITH COMMA: (e.g., 1,5,9,10,14,15)

The criteria for selecting a peak are: (1) no significant artifact should be present in this waveform (and in the waveform immediately preceeding it), (2) delta-t must be correctly marked, and (3) waveforms must occur during the exhalation respiratory plateau (pause between breaths). Referring to figure 1, the user would likely select peaks 1,5,9,10,14, and 15.

After making peak selections, the user presses the "RETURN" key and is prompted with the query

DO YOU WANT A COPY [Y/N]:

Answering "Y" will produce a hard copy of the epoch being displayed on the screen. Answering "N" will prompt the query

DO YOU WANT TO EXIT [Y/N]:

Entering "N" causes the program to display the next 15-sec epoch of data from which the user may select peaks for data reduction. This is repeated until all 15-sec data epochs of this file have been processed. Note that if the user has just completed the last epoch of data in this file, the program terminates automatically (i.e., responds as though a "Y" was entered). Entering "Y" terminates the data reduction phase of the program and the screen will display options for user selection:

ORIGINAL PAGE IS
OF POOR QUALITY

ENTER 1 TO SAVE DATA ON FLOPPY AND EXIT

ENTER 2 TO EXIT (SAVE DATA ON HARD DISK)

ENTER 3 TO RESTART DATA EDITNG

ENTER 4 TO RESTART CALIBRATION ACQUISITION

[ENTER 1, 2, 3 OR 4]:

Entering "1" will cause the program to prompt the user for additional key field information used to identify this file. The screen displays

GROUP TYPE (maximum of 8 characters): (e.g., CONTROL)

RUN NUMBER: (e.g., 1)

DIRECTION (cc,cw,nd): NOTE: cc=counterclockwise, cw=clockwise
nd=no direction (i.e., no rotation)

SUBJECT'S INITIALS (first and last): (e.g., PC)

AGE: (e.g., 49)

SUSCEPTIBILITY (1,2,3): NOTE: 1 = high motion sickness
susceptible, 2 = moderate and
3 = low.

SEX (m/f): (e.g., M)

TEST DATE (mmddyy): (e.g., MAY0288)

TEST TIME (military in hr, min, e.g., for 1:00 pm enter 1300):

PRE HEART RATE (in beats/min): (e.g., 68)

PRE TEMPERATURE: (e.g., 97.6)

PRE B.P. (sys,dia, e.g., 120,80):

PRE/POST TEST BASELINE (minutes): (e.g., 10)

PLEASE INSERT DATA FLOPPY WITH AT LEAST 60 CONTIGUOUS FREE BLOCKS
INTO DRIVES DZ1 & DZ2. WHEN READY, PRESS RETURN

ARE YOU FINISHED WITH DIGITIZING DATA? [Y/N]

Entering "Y" will prompt the user as follows:

DO YOU WANT TO DELETE THE DIGITIZED FILE? [Y/N]:

(NOTE: After this input, the program returns
the user to the original menu).

Entering "2" will return the user to the original menu for this
program.

Entering "3" will return the user to the editing portion of
this program.

Entering "4" will return the user to the calibration acquisition
portion of this program.

REFERENCES

1. Cowings, P.S.; Suter, S.; Toscano, W.B.; Kamiya, J.; and Naifeh, K.: General Autonomic Components of Motion Sickness, Psychophysiology, vol. 23, no. 3, 1986, pp. 542-551.
2. McGregor, M.; Donevan, R.E.; and Anderson, N.M.: Influence of Carbon Dioxide and Hyperventilation on Cardiac Output in Man. J. Appl. Physiol., vol. 17, no. 6, 1962, pp. 933-937.
3. Miller, J.C.; and Horvath, S.M.: Impedance Cardiography, Psychophysiology, vol. 15, no. 5, 1978, pp. 80-91.

APPENDIX

SOURCE CODE

```

C      PROGRAM IMPMLT
C
C      AN IMPEDANCE CARDIOGRAPH PROGRAM USING DIGITIZED DATA FILES
C
C      PARAMETER (IFREQ=200)
C      DIMENSION SLOPE(2),ENTRCP(2),CALVAL(2,2)
C      COMMON ISPACE(4*IFREQ ),SPACE(30*IFREQ )
C      REAL*8      DRAW ,IMPACQ ,IMPCAL,IMPSTR
C
C      =====
C      DATA ACQUISITION OF BIOFEEDBACK PARAMETERS
C      =====
C
C      PEAL TEMPST,
C      @          ZERO,CFREQ,TFREQ,EPOCH
C
C      INTEGER PREDIA,PREHR,PRETEM,PRESYS,PSTDIA,PSTSYS,NUMDIA,
C      @          PSTHR,PSTTEM,RUNNUM,IDAY,NSN,AGE,SUSC,
C      @          TESTIM,MIN,PSTTST,PART,CHANEL,TOTMIN,PRETST,
C      @          IMINS,IAN,          HIGH( 2,2)          ,DIAG(15),
C      @          CLEANS(2),ICLEAN(2)
C
C      LOGICAL*1 SEX,ANS,TEMP,DMY,IRSP,IABORT,IPAWS,TYMEAN,STDDEV,
C      @          GROUP(8),DATE(7),DIRECT(2),OSN(2),STRNG(8),
C      @          FRMT(12)
C      CHARACTER*1 CMT(80),RUN(4),RUNTYP,DATFIL(14),EXP
C      CHARACTER*14 CATFIL
C      EQUIVALENCE (DATFIL,CATFIL)
C
C
C      -----
C      INITIALIZATION
C      -----
C
C      DATA IMPCAL/'IMPCAL'/ ,IMPSTR/'IMPSTR'/
C      DATA DRAW/'DRAWRM'/ ,ICLEAN/27,99/,IMPACQ/'IMPACQ'/
C      DATA CLEANS/27,99      /,CATFIL/'00000000000000'/
C
C      TYPE 500,CLEANS
500    FORMAT(X,4A1)
C      TYPE 101
101    FORMAT(10X,60('*'),//,27X,' IMPEDANCE CARDIOGRAPH PROGRAM',
1    //,10X,60('*'),///// )
125    CONTINUE
C      CATFIL(1:14)='00000000000000'
C      TYPE 102

```

```

102     FORMAT(20X,45('*'),/,20X,'*',
1      /,20X,'*',6X,' SELECT FROM THE FOLLOWING',/,20X,'*',
1      /,20X,'*',10X,'1. DIGITIZE A DATA FILE FROM TAPE',
2      /,20X,'*',10X,'2. DATA REDUCTION OF DIGITIZED FILE',
3      /,20X,'*',10X,'3. EXIT'.
4      /,20X,45('*'),/,
5      /' [ENTER 1, 2 OR 3]: '$)
150    CONTINUE
        ACCEPT * , IANS
        IF(IANS.EQ.1) GOTO 200
        IF(IANS.EQ.2) GOTO 2000
        IF(IANS.EQ.3) GOTO 99999
        TYPE *, ' ***** E R R O R - TYPE 1 OR 2 OR 3 ONLY. NOW DO IT'
        GOTO 150
200    CONTINUE
        IEFN=6
C      CALL IMPACQ(CATFIL)
        CALL SPAWN(RAD50(IMPACQ),,,IEFN,,IESD,,,,,IDS)
D      TYPE *, ' SPAWN CALLED ',IEFN= ',IEFN,' IESD= ',IESD,'IDS= ',IDS
        CALL WAITFR(IEFN,IDS2)
C      TYPE *, ' IEFN SET NOW ',IDS2
        TYPE 500,CLEANS
        GOTO 125
2000   CONTINUE

C      =====
C      PHASE 1 -- KEY FIELD ENTRY
C      =====

        TYPE 500,CLEANS
        TYPE 501
501    FORMAT(10X,60('*'),///// ,27X,' PHASE 1 -- KEY FIELD ENTRY',///// ,
1      10X,60('*'),///// )

        TYPE 3998
3998   FORMAT (/ ' ENTER RUN ID FOR THIS FILE:',$)
        ACCEPT 3999,CATFIL(5:10)
3999   FORMAT(A6)
        TYPE 4000
4000   FORMAT(/ ' LENGTH (CM.) BETWEEN INNER ELECTRODES: ', $)
        ACCEPT * ,ELEN
        TYPE 4100
4100   FORMAT(/ ' HERMATOCRIT COUNT : ', $)
        ACCEPT * ,HEM
        TYPE 4200
4200   FORMAT(/ ' HEIGHT(CM) : ', $)
        ACCEPT * ,HEIGHT
        TYPE 4300
4300   FORMAT(/ ' WEIGHT(KG) : ', $)
        ACCEPT * ,WEIGHT

```

```

C      -----
C      SET DEFAULTS AND ZERO CURRENTLY UNUSED FIELDS
C      -----

```

ORIGINAL PAGE IS
OF POOR QUALITY

```

      RHO=53.2*EXP      (.022*HEM)
      BSA=HEIGHT** .725*WEIGHT** .425*.007184
      C=RHO*ELEN**2
D     TYPE *, ' C= ', C, ' RHO= ', RHO
D     CALL WAIT(5,2)

C     =====
C     PHASE 2 -- CALIBRATION VALUES
C     =====

6000  CONTINUE
      CALL ERRSET(59,.TRUE.,.FALSE.,.TRUE.,.FALSE.,MAX)
      CALL IMPVAL (CALVAL)

C     =====
C     PHASE 3 -- CALIBRATION ACQUISITION
C     =====

7010  CONTINUE
      OPEN (UNIT=1,NAME='DW1:NUDRAW.TEL',TYPE='NEW',FORM='FORMATTED')
      WRITE(1,2) CALVAL      ,C,BSA,CATFIL,IMINS
2     FORMAT(6( G15.7),A14,I10)
      CLOSE (UNIT=1)
      IEFN=6
      CALL SPAWN (RAD50(IMPVAL),,,IEFN,,IESD,,,,,IDS)
D     TYPE *, ' SPAWN CALLED ', 'IEFN= ',IEFN, ' IESD= ',IESD, 'IDS= ',IDS
      CALL WAITFR(IEFN,IDS2)
C     TYPE *, ' IEFN SET NOW ',IDS2
      OPEN (UNIT=1,NAME='DW1:NUDRAW.TEL',TYPE='OLD',FORM='FORMATTED')
      READ (1,2) SLOPE,ENTRCP,C,BSA,CATFIL,IMINS
      CLOSE (UNIT=1)

      CATFIL(11:14)='.IMP'

C
C     OPEN INTERNAL FILE THAT STORES RESULTS OF THE REDUCTION PROCESS
C     IF IT EXISTS AND FIND OUT WHAT EPOCHS HAVE BEEN COMPLETED. OTHER
C     WISE CREATE A NEW INTERNAL FILE TO BEGIN REDUCTION
C

      OPEN (UNIT=1,STATUS='OLD',NAME=CATFIL,
1     FORM='UNFORMATTED',ACCESS='DIRECT',RECL=8,ERR=50)
      CLOSE(UNIT=1)
      CALL INVEST(CATFIL) !DETERMINE WHAT EPOCHS HAVE BEEN COMPLETED
      TYPE *, 'PRESS [RETURN] TO CONTINUE '
      ACCEPT 500 ,ANS
      GOTO 55
50  CONTINUE
      OPEN (UNIT=1,STATUS='NEW',NAME=CATFIL,
1     FORM='UNFORMATTED',ACCESS='DIRECT',RECL=8)
      CLOSE(UNIT=1)
      MAXREC=380
      IF(CATFIL(10:10).EQ.'0') MAXREC=120
      IF(CATFIL(10:10).EQ.'4') MAXREC=220
      CALL ZEREC (MAXREC,CATFIL) !ZERO CONTENTS OF NEW FILE FOR MAX RECS
55  CONTINUE

```

ORIGINAL PAGE IS
OF POOR QUALITY

```

C      IEFN=6
C
C      =====
C      PHASE 4 -- DATA REDUCTION
C      =====
C
C      NOTE FOR THIS SPAWNING TO WORK TYPE INSTALL DRAWRM.TSK TO
C      INSTALL THE DRAWED.FTN PROGRAM
C
80000  CONTINUE
      CALL SPAWN(RAD50(DRAW),,,,IEFN,,IESD,,,,,IDS)
D      TYPE *, ' SPAWN CALLED ', 'IEFN= ', IEFN, ' IESD= ', IESD, 'IDS= ', IDS
      CALL WAITR(IEFN,IDS2)
C      TYPE *, ' IEFN SET NOW ',IDS2
C      =====
C      PHASE 5 -- DATA STORAGE
C      =====
C
C
9000  CONTINUE
      TYPE 500,ICLEAN

9101  CONTINUE
      TYPE 9109
9109  FORMAT(////////////////)
      TYPE 9110
9110  FORMAT(' WOULD YOU CARE TO ',
1      /7X,'1) SAVE DATA ON FLOPPY AND EXIT',
2      /7X,'2) EXIT (SAVE DATA ON DISK)',
3      /7X,'3) RESTART EDITING',
4      /7X,'4) RESTART CALIBRATION ACQUISITION (SAVE)',
5      /' [ENTER 1,2,3,OR 4 ]: ', $)
      ACCEPT 10,TEMP
10    FORMAT(A1)
      IF(TEMP.EQ.'2'.OR .TEMP.EQ.'4') THEN
      OPEN(UNIT=1, NAME='DW1:NUDRAW.TEL',STATUS='OLD')
      CLOSE (UNIT=1,DISPOSE='DELETE')
      END IF
      IF (TEMP.EQ.'1'.OR.TEMP.EQ.'2'.OR.TEMP.EQ.'3'.OR.TEMP.EQ.'4')
1      GO TO 9114
      TYPE 9108
9108  FORMAT(///// ' YOU MUST NOT HAVE CHOSEN A VALID ALTERNATIVE;',
@      ' CHOOSE ONE OF THE FOLLOWING: '/')
      GOTO 9101
9114  CONTINUE
      DECODE (1,110,TEMP) ITEMP
110   FORMAT(I1)
      GOTO (9115,99990,80000,6000) ITEMP

```


OF POOR QUALITY

```

C -----
C ASSEMBLE DATAFILE NAME
C -----

9115    CONTINUE
        CALL SPAWN(RAD50(IMPSTR),,,,IEFN,,IESD,,,,,IDS)
D      TYPE *,' SPAWN CALLED ','IEFN= ','IEFN,' IESD= ','IESD,'IDS= ','IDS
        CALL WAITFR(IEFN,IDS2)
C      TYPE *,' IEFN SET NOW ',IDS2
C      CALL IMPSTR(CATFIL,IFLAG,GROUP,RUNNUM,DIRECT,OSN,AGE,SUSC,SEX,
C 1     DATE,TESTIM,PREHR,TEMPRT,PRETEM,PRESKS,PREDIA,PRETST,PSTTST)
C      TYPE 500,CLEANS
        OPEN(UNIT=1, NAME='DW1:NUDRAW.TEL',STATUS='OLD')
        CLOSE (UNIT=1,DISPOSE='DELETE')
        GOTO 125
C10000  CONTINUE
C      IFLAG=1 ! WRITE COMMENTARY FILE
C      CALL COMENT (KCHAR)
C      GOTO 9101
99990   CONTINUE
        TYPE 500,CLEANS
        GOTO 125
99999   CONTINUE
        END
        SUBROUTINE ZEREC(MAXREC,CATFIL)
        integer tcount
        CHARACTER*14 CATFIL
        OPEN (UNIT=1,STATUS='OLD',NAME=CATFIL,
1       FORM='UNFORMATTED',ACCESS='DIRECT',RECL=8)

        do 100 tcount = 1,MAXREC
            write(1,REC=TCOUNT,ERR=999)0.,0.,0.,0.,0.,0.,0.,0.
100      continue

999     CLOSE(UNIT=1, STATUS='SAVE')
        RETURN
        END
        SUBROUTINE INVEST(CATFIL)
        DIMENSION BUF(8,1),NONZER(400)
        CHARACTER*14 CATFIL
        CALL ERRSET(36,.TRUE.,.FALSE.,.TRUE.,.FALSE.)
        IREC=1
        INUM=0
        TYPE *,' FROM THE RECORD, THE FOLLOWING EPOCHS ARE NONZERO :'
100     CONTINUE
        OPEN (UNIT=1, status='old', NAME=CATFIL,ERR=999,
1       FORM='UNFORMATTED', access='direct',recl=8 )

            read(1,REC=IREC ,ERR=999)(BUF(ICOUNT,1),ICOUNT=1,8)
            IF( BUF(1,1)) 999,999,200
200     CONTINUE
            INUM=INUM+1
            NONZER(INUM)=IREC
999     CLOSE(UNIT=1, STATUS='SAVE')

```

```

CALL ERRNST(36,IERR)
IF(IERR.EQ.1) THEN
TYPE *,(NONZER(1),I=1,INUM)
RETURN
ELSE
IREC=IREC+1
GOTO 100
END IF
END
PROGRAM IMPCAL
C (SLOPE,INTRCP,CATFIL)
PARAMETER (IFREQ=200)
REAL CALVAL(2, 2),INTRCP(2),SLOPE(2)

INTEGER CLEANS(2),CHANEL,HIGH(2,2)
COMMON ISPACE(4*IFREQ),SPACE(30*IFREQ)
EQUIVALENCE(CALVAL,SPACE(1))
CHARACTER*14 CATFIL
DATA CLEANS/27,99/
C =====
C PHASE 3 -- CALIBRATION ACQUISITION
C =====
CFREQ=IFREQ
OPEN(UNIT=1,NAME='DW1:NUDRAW.TEL',TYPE='OLD',ORM='FORMATTED')
READ (1,2)CALVAL,C,BSA,CATFIL,IMINS
2 FORMAT(6(G15.7),A14,I10)
CLOSE(UNIT=1)
7010 TYPE 500,CLEANS
500 FORMAT(X,4A1)
TYPE 700
700 FORMAT(10X,60('*'),/////,22X,' PHASE 3 -- CALIBRATION
1 ACQUISITION',/////,10X,60('*'),/////)
CATFIL(1:4)='DW1:'
CATFIL(11:14)='.DIG'

C -----
C CALIBRATE BASELINE IMPEDANCE - CHANNEL 2
C -----
CALL CALZ (SLOPE(1),INTRCP(1),CALVAL,CATFIL)

C -----
C CALIBRATE dZ/dT SIGNAL - CHANNEL 3
C -----
CALL CALDZT (SLOPE(2),INTRCP(2),CALVAL,CATFIL)

OPEN(UNIT=1,NAME='DW1:NUDRAW.TEL',TYPE='OLD',FORM='FORMATTED')
WRITE(1,2)SLOPE,INTRCP,C,BSA,CATFIL,IMINS
CLOSE(UNIT=1)
END
SUBROUTINE DIRECT (TCOUNT,CATFIL,NCHAN)
C STORAGE SEQUENCE AS FOLLOWS:
C FIRST - RAW RESPIRATION

```

ORIGINAL PAGE IS
OF POOR QUALITY

OF POOR QUALITY

C SECOND - DELTA Z
C THIRD - DZ/DT
C FOURTH - ECG

c TKB FOR TASK BUILD

C
c TKB>FILENAME=FILENAME
C TKB>/
C TKB>MAXBUF=1600
C TKB>//
c

PARAMETER (IFREQ=200)
integer buf(IFREQ), tcount
DIMENSION X(15*IFREQ),Y(15*IFREQ)
COMMON BUF,X,Y
CHARACTER*14 CATFIL

1 OPEN (UNIT=1, status='old', NAME=CATFIL
FORM='UNFORMATTED', access='direct',recl=4*IFREQ)

100 read(1,REC=TCOUNT,ERR=999)(BUF(ICOUNT),ICOUNT=1,4*IFREQ)
continue

999 CLOSE(UNIT=1, STATUS='SAVE')
10 FORMAT(x,i3,4(2X,i4,2x,i4))
20 FORMAT(1X,i5)

RETURN

END

SUBROUTINE CALZ (SLOPE,INTRCP,CALVAL,CATFIL)

C -----
C CALIBRATE BASELINE IMPEDANCE - CHANNEL 2
C -----
C

REAL CALVAL(2, 2),CFREQ,TFREQ,EPOCH,INTRCP
INTEGER*4 IADSUM
INTEGER CHANEL,ICOUNT,HIGH(2),
@ CLEANS(2),MAXPTS,ISTAT(2),IEFN,MDSYN,
@ ICHAN(8),ICONV,IFORM,ITRIG,ISTAT2(2),IEFN2,NPTS(2)

CHARACTER*4,IWORD(2)
CHARACTER*14,CATFIL

PARAMETER (IFREQ=200)
COMMON ISPACE(4*IFREQ),SPACE(30*IFREQ)
LOGICAL*1 ANS

DATA IWORD /'HIGH','LOW '/
DATA CLEANS/27,99/

C
C
C

```

C      -----
C      INITIALIZATION
C      -----
C      =====
C      PERFORM PRELIMINARY ACQUISITION
C      =====
C      CALL WAIT( 2,2)
C      TYPE 500,CLEANS
500    FORMAT(X,4A1)
C      TYPE 700
700    FORMAT(10X,60('*'),/////.22X,' CALIBRATION
;      OF BASELINE IMPEDANCE',/////,10X,60('*')////////)

C      CALL WAIT( 2,2)
C      CHANEL=2
C      -----
C      -----
7001   CONTINUE
C      -----
C      -----
C      -----
C      OPEN(UNIT=2,NAME='SY:DUMP.TST',FORM='FORMATTED',TYPE='NEW')
C      ISTREC=1
7019   CONTINUE
C      N=1 !HIGH CALIBRATIONS
C      -----
C      PERFORM HIGH CALIBRATION
C      -----
C      CALIB=.TRUE.
7110   CONTINUE
C      -----
C      INITIALIZE SLOPE, INTERCEPT
C      -----

C      SLOPE = 0.0
C      INTRCP = 0.0
C      NSAMP=0
C      IADSUM=0

C      DO 7000 I=0,29
C      IREC= 30*(N-1)+ISTREC+I
C      CALL DIRECT(IREC,CATFIL)
C
C
C
C      -----
C      CALCULATE MEAN
C      -----

```

ORIGINAL PAGE IS
OF POOR QUALITY

```

NPTCHN=IFREQ !NUMBER OF SAMPLES PER CHANNEL PER RECORD=PER SECOND
DO 704 J=1,NPTCHN
  IADSUM=IADSUM+ISPACE ( (J-1)*4+CHANNEL )
C   TYPE *,IADSUM,ISPACE ( (J-1)*4+CHANNEL )
  NSAMP=NSAMP+1
C   WRITE(2,19999)IAD(I+(J-1)*CHANNEL+(K-1)*MAXPTS),NPTCHN,CHANNEL,
C   1 NPTS(K),NSAMP,I,J,K,N
C19999  FORMAT(9I8)
704    CONTINUE
C   WRITE(2,*) IADSUM,NSAMP,ISPACE
7000   CONTINUE
      HIGH(N)=IADSUM/NSAMP

C
C   -----
C   ECHO MEAN AND VERIFY
C   -----
C

      TYPE 715,IWORD(N)
715   FORMAT(////////' CHANNEL ',5X' SIGNAL ',5X,'A/D',A4,
1      /,' -----',5X,'-----',5X,'-----')

      TYPE 71550 , HIGH(N)
71550  FORMAT( 5X,' 1',8X,'IMPEDANCE ',5X,I10)
      TYPE 717,IWORD(N)
717    FORMAT(/' DO YOU WANT TO RE-RUN ',A4,' CALS [Y/N]? ', $)
      ACCEPT 617,ANS
617    FORMAT(A1)
      IF (ANS .EQ. 'Y'.OR.ANS.EQ.'y') THEN
        TYPE *, ' WHAT STARTING RECORD '
        ACCEPT *,ISTREC
        GOTO 7110
      END IF
      IF (N.EQ.2) GOTO 724
      N=2
      GOTO 7110

724    CONTINUE

C   -----
C   CALCULATE SLOPE AND INTERCEPT
C   -----

      IF (HIGH(1) .NE. HIGH(2)) GOTO 7185
      INTRCP = 0.0
      SLOPE = 0.0
      GOTO 718
7185   SLOPE = (CALVAL(2,1)-CALVAL(1,1))/FLOAT (HIGH(1)-HIGH(2))
      INTRCP = CALVAL(1,1) - SLOPE*FLOAT(HIGH(2))
718    CONTINUE

C   -----
C   ECHO SLOPE,INTERCEPT
C   -----

```

```

      TYPE 720
720  FORMAT(/' CHANNEL  ',5X,' SIGNAL  ',5X,' ADLOW  ',5X,
      1' ADHI    ',5X,
      2 ' LOWCAL ',5X,' HICAL  ',/,X,6('-----',5X))
      TYPE 81550 , HIGH(2)
81550  FORMAT( 4X,' 1',7X,'IMPEDANCE ',4X,I10,$)
      TYPE 722,HIGH (1),CALVAL(1,1),CALVAL (2,1)
721  CONTINUE
722  FORMAT('+',4X,I10,4X,F10.4,4X,F10.4)
      TYPE 723
723  FORMAT(/' WOULD YOU LIKE TO REPEAT CALIBRATIONS [Y/N]?',$)
      ACCEPT 617,ANS
      IF (ANS .EQ. 'Y'.OR.ANS.EQ.'y') THEN
      TYPE *, ' WHAT STARING RECORD '
      ACCEPT *,ISTREC
      GOTO 7019
      END IF
      CLOSE (UNIT=2)
      RETURN
      END
      FUNCTION MAXIAD (IAD,MAXPTS,NPTS)
      DIMENSION IAD(NPTS*MAXPTS)
      MAXIAD= IAD(3)
      DO 100 J=1,MAXPTS
      IF(IAD((J-1)*NPTS+3 ).LT.MAXIAD) GOTO 100
      MAXIAD= IAD((J-1)*NPTS+3)
100  CONTINUE
      RETURN
      END
      FUNCTION MINIAD (IAD,MAXPTS,NPTS)
      DIMENSION IAD(NPTS*MAXPTS)
      MINIAD= IAD(3)

      DO 100 J=1,MAXPTS
      IF(IAD((J-1)*NPTS+3 ).GT.MINIAD) GOTO 100
      MINIAD=IAD((J-1)*NPTS+3)
100  CONTINUE
      RETURN
      END
      SUBROUTINE CALDZT (SLOPE,INTRCP,CALVAL,CATFIL)
      C -----
      C CALIBRATE dZ/dT SIGNAL - CHANNEL 3
      C -----
      C
      REAL CALVAL(2, 2),CFREQ,TFREQ,EPOCH,INTRCP
      INTEGER*4 IADSUM,IAD4,IAD4MX,IADELT
      INTEGER CHANEL,ICOUNT,HIGH(2),IAD(400),
      @      CLEANS(2),MAXPTS,ISTAT(2),IEFN,MDSYN,
      @      ICHAN(8),ICONV,IFORM,ITRIG,ISTAT2(2),IEFN2,NPTS(2)

      PARAMETER (IFREQ=200)
      COMMON ISPACE(4*IFREQ),SPACE(30*IFREQ)
      EQUIVALENCE (IAD,ISPACE)

```

```

CHARACTER*4,IWORD(2)
CHARACTER*14,CATFIL

C
C
C
C
C  -----
C  INITIALIZATION
C  -----

DATA CLEANS/27,99/
C  =====
C  PERFORM PRELIMINARY ACQUISITION
C  =====
TYPE 500,CLEANS
500  FORMAT(X,2A1)
TYPE 700
700  FORMAT(10X,60('*'),/////,22X,' CALIBRATION
1    OF dZ/dT SIGNAL',/////,10X,60('*'),/////)
CALL WAIT(5,2)

CHANEL=3
C  -----
C  -----

7001  CONTINUE
C  -----
C  -----
C
C
7015  CONTINUE
C  -----
C  -----
C
C
7019  CONTINUE
C  OPEN(UNIT=2,NAME='SY:DUMP.TST',FORM='FORMATTED',TYPE='NEW')

N=1 !HIGH CALIBRATIONS
SLOPE=0.0
INTRCP=0.0
713  CONTINUE
NCHAN=0
NSAMP=0
IADSUM=0

C  -----
C  READ IN A-D VALUES
C  -----

IREC= 31+I
C  TYPE *, IREC= ',IREC,'N= ',N,'I= ',I
CALL DIRECT(IREC,CATFIL)
D  TYPE *,(IAD((J-1)* 4 +3),J=1,100)

```

ORIGINAL PAGE IS
OF POOR QUALITY

```

73100 CONTINUE
71300 CONTINUE
C
C
C
C
C -----
C CALCULATE MEAN
C -----
C
C IF(N.EQ.1) THEN
C FIND MAX VALUE IF N=1
C IADMAX=MAXIAD(IAD,IFREQ,4 )
D TYPE *, ' IADMAX= ', IADMAX
C IADELT=IADMAX
C IADELT=IADELT*5
C IF(IADELT.LT.0) IADELT=-IADELT
C ELSE
C FOR N=2 FIND MAX AND MIN
C IADMAX=NUL(IAD,MAXPTS,NPTS)
C IADMAX=MAXIAD(IAD,IFREQ,4 )
C IADMIN=MINIAD(IAD,IFREQ,4 )
C IZERO=IADMIN+(IADMAX-IADMIN)*2/5
C END IF
C NPTCHN=IFREQ
C DO 704 J=1,NPTCHN
C IF(N.EQ.1) THEN
C USE IAD VALUE IF IAD >= MAX VALUE- 5%*(MAX VALUE)
C
C IAD4=IAD((J-1)*4+CHANEL )
C IAD4MX=IADMAX
C IF((100*IAD4).GE.(IAD4MX*100-IADELT )) THEN
C NCHAN=NCHAN+1
C IADSUM=IADSUM+IAD ((J-1)*4+CHANEL )
C TYPE *,IAD4,IADELT,IAD4MX
C END IF
C ELSE
C IF(IAD((J-1)*4+CHANEL ).LE.IZERO +80.AND.
1 IAD((J-1)*4+CHANEL ).GE.IZERO -80) THEN
C NCHAN=NCHAN+1
C IADSUM=IADSUM+IAD ((J -1)*4+CHANEL )
C END IF
C END IF
C WRITE(2,19999)IAD(J+(K-1)*MAXPTS),NPTCHN,CHANEL,
C 1 NPTS(K),NCHAN,I,J,K,IADMAX
19999 FORMAT(9I8)
704 CONTINUE
7000 CONTINUE
C HIGH(N)=IADSUM/NCHAN
C IF (N.EQ.2) GOTO 714
C N=2
C GOTO 713
714 CONTINUE

```



```

C      -----
C      ECHO MEANS AND VERIFY
C      -----
C
C      WRITE(2,*      ) IAD,HIGH,IADMAX,IADMIN,IZERO
C      TYPE 715,IWORD(N)
715   FORMAT(////////' CHANNEL  ',5X'   SIGNAL  ',5X,'A/D',A4,
1      /,' -----',5X,'-----',5X,'-----')
C      TYPE 71551 , HIGH(N)

71551  FORMAT( 5X,' 2',8X,'   dZ/dT ',5X,I10)

C      TYPE 717
717   FORMAT('/ DO YOU WANT TO RE-RUN dZ/dT CALS [Y/N]? ',,$)
C      ACCEPT 617,ANS
617   FORMAT(A1)
C      IF (ANS .EQ. 'Y'.OR.ANS.EQ.'y') GOTO 7110

      CLOSE(UNIT=2)

C      -----
C      CALCULATE SLOPE AND INTERCEPT
C      -----

      IF (HIGH(1) .NE. HIGH(2)) GOTO 7135
      INTRCP = 0.0
      SLOPE = 0.0
      GOTO 718
7185  SLOPE = (CALVAL(2,2)-CALVAL(1,2))/FLOAT (HIGH(1)-HIGH(2))
      INTRCP = CALVAL(1,2) - SLOPE*FLOAT(HIGH(2))
718   CONTINUE

C      -----
C      ECHO SLOPE, INTERCEPT
C      -----

      TYPE 720
720   FORMAT('/ CHANNEL  ',5X,' SIGNAL  ',5X,' ADLOW  ',5X,
1' ADHI    ',5X,
2' LOWCAL ',5X,' HICAL  ',/,X,6('-----',5X))
      TYPE 81551 , HIGH(2)
81551  FORMAT( 4X,' 2',7X,'   dZ/dT ',4X,I10,$)
      TYPE 722,HIGH (1),CALVAL(1,2),CALVAL (2,2)
722   FORMAT('+',4X,I10,4X,F10.4,4X,F10.4)
      TYPE *, ' SLOPE ',SLOPE,' INTERCEPT ', INTRCP
      TYPE 723
723   FORMAT('/ WOULD YOU LIKE TO REPEAT CALIBRATIONS [Y/N]? ',,$)
      ACCEPT 617,ANS
      IF (ANS .EQ. 'Y'.OR.ANS.EQ.'y') THEN
      TYPE *, ' WHAT STARING RECORD '
      ACCEPT *,ISTREC
      GOTO 7019
      END IF

```



```

C      =====
C      PHASE 5 -- DATA ACQUISITION
C      =====
      FREQ=IFREQ
      MAXPTS=4*IFREQ
      CALL GETADR(IPAR(1),IBUF)
      IPAR(2)=1
9000   TYPE 500,CLEANS
      TYPE 900
900    FORMAT(10X,60('*'),///// ,22X,' DIGITIZING A DATA FILE ',
1      ///// ,10X,60('*'),///// )

C      -----
C      REQUEST MAXIMUM RUN DURATION
C      -----

      TYPE 90109
90109  FORMAT(///' PLEASE ENTER RUN ID FOR THIS DATA FILE: ', $)
      ACCEPT 618,(CATFIL(5:10))
      CATFIL(1:4)='DW1:'
      CATFIL(11:14)='.DIG'
      TYPE 90110
90110  FORMAT(///' PLEASE ENTER RUN TIME (IN SECS.)FOR TAPE FILE: ', $)
      ACCEPT *,MAXRUN
      ISIZE=6.25 *MAXRUN+1
      TYPE 90111
90111  FORMAT(///' HOW FAST TO SAMPLE THE DATA ON TAPE ? ',
1      /10X,' ENTER 1 FOR REAL TIME RATE ',
2      /10X,' ENTER 2 FOR TWICE REAL TIME ',
3      /' [ENTER 1 OR 2] NOW: ', $)
      ACCEPT *,ISAMP

C      NOTE: MAXPTS=FREQ*4
      CFREQ=FREQ*ISAMP
      CALL SETFRQ(CFREQ,TFREQ)
      DO 190 I=1,4
      ICHAN(I)=1
190    CONTINUE
      IEFN1=8
      IEFN2=10
      ICOUNT=2*MAXPTS

C      -----
C      OPEN TEMPORARY DATA FILE
C      -----

9015   CONTINUE
C      OPEN(UNIT=2,NAME='SY:TEMP.DAT',TYPE='NEW')
      OPEN (UNIT=1,STATUS='NEW',NAME=CATFIL,INITIALSIZE=-ISIZE,
1      FORM='UNFORMATTED',ACCESS='DIRECT',RECL=4*IFREQ)

      MIN = 0
      IMINS = 0

```

```

C      -----
C      INITIALIZE SUM,ICOUNT
C      -----

9011  CONTINUE
      ICONT(1)=ICOUNT/2
      ICONT(2)=ICOUNT/2
      CALL AINIT (ISTAT,IEFN1)
      IF(ISTAT(1).EQ."40000.OR.ISTAT(1).EQ."40040) GOTO 7000
      TYPE *, ' ISTAT(1) = ',ISTAT(1)
      PAUSE
7000  CONTINUE

      TYPE *, ' '
      IF (IRSP .NE. IPAWS)GO TO 9022
      PAUSE
      TYPE 500,BLINK
      TYPE 7100
7100  FORMAT(/' *****',
1      /' * DATA COLLECTION RESUMED *',
2      /' *****')
      TYPE 500,RESET
      GO TO 901
9022  TYPE 9012
9012  @  FORMAT (' HIT "S" TO START, "P" TO PAUSE,',
      ' "A" TO ABORT: '/' FOLLOWED BY A [RETURN]')
      ACCEPT 617, ANS
      IF (ANS .EQ. 'A') GO TO 9100
      IF (ANS .EQ. 'P') PAUSE
      TYPE 500,CLEANS
      TYPE 500,BLINK
      TYPE 97100
97100 FORMAT(/////////,
1      /' *****',
2      /' * DATA COLLECTION UNDERWAY *',
3      /' *****')

      TYPE 500,RESET
901  CONTINUE
      CALL QIO(IREAD,7, 4,,IOSB,IPAR,IDSW)
      CALL AIN (ISTAT, ! START CLOCK
2  IAD ,
3  ICONT(1),
4  IEFN1,
5  MDSYN,
6  ICHAN,
7  ICONV,
8  IFORM,
9  ITRIG,ITIME,CFREQ)
      CALL ABUF(ISTAT2,IAD2,ICONT(2),IEFN2,MDDE ,ITIME)
908  CONTINUE
      IMINS = IMINS + 1
C      TYPE *, ' IMINS= ',IMINS
      K=1

```

```

C      N=1
C      IEND=ICOUNT
80001  CONTINUE
      IISTAT=ISTA(1+2*(K-1))
      IF (IISTAT .EQ."40000.OR.IISTAT .EQ."40040)
1      GOTO 81300
C      IF (IISTAT.NE.0) TYPE *,K,IISTAT
      TYPE *,K,IISTAT,IMINS
      GOTO 80001
81300  CONTINUE
C
C
C      -----
C      STORE DIGITIZED DATA FOR EACH SECOND
C      -----
C
      CALL DIRECT(IMINS,K)
      NCONT=ICONT(K)
      CALL ABUF(ISTA(1+2*(K-1)),IAD(1+(K-1)*MAXPTS),NCONT ,IEFN(K),
1  MODE,ITIME)

C      -----
C      ECHO MEANS AND STANDARD DEVIATIONS
C      -----

C      -----
C      IF END OF SESSION, END DATA ACQUISITION
C      -----
      MIN = IMINS
      IF(MIN.EQ. 30) THEN
      CALL ACSTAT(ISTAT,JSTAT,1,IEFN1,0)
      CALL QIO("12,9)
      CALL QIO("12,7)
      PAUSE ' - ADVANCE TAPE TO START OF LOW CALS NOW; THEN '
      CALL AINIT(ISTAT,IEFN1,0)
      TYPE 500,BLINK
      TYPE 7100
      TYPE 500,RESET
      GOTO 901
      END IF
      IF(MIN.EQ. 60) THEN
      CALL ACSTAT(ISTAT,JSTAT,1,IEFN1,0)
      CALL QIO("12,9)
      CALL QIO("12,7)
      PAUSE ' - ADVANCE TAPE TO START OF DATA NOW; THEN '
      CALL AINIT(ISTAT,IEFN1,0)
      TYPE 500,BLINK
      TYPE 7100
      TYPE 500,RESET
      GOTO 901
      END IF
      IF (MIN .GE. MAXRUN) GO TO 9100

```

```

C      -----
C      IF PAUSE OR ABORT, END DATA ACQUISITION
C      -----
      CALL READEF( 4,IEF)
      IF(IEF .GT.0) THEN
      IRSP =IBUF
      IF (IRSP .EQ. IPAWS) THEN
          CALL ACSTAT(ISTAT,JSTAT,1,IEFN1,0)
      CALL QIO("12,9)
C      CALL QIO("12,7)
      CALL AINIT (ISTAT,IEFN1,0)
      GOTO 7000
      ELSE
      IF (IRSP .EQ. IABORT) GO TO 9100
      END IF
C      SOMETHING ELSE TYPED IN; QUE I/O REQUEST AGAIN
      CALL QIO(IREAD,7, 4,,IOSB,IPAR,IDSW)
      END IF
      IK=K
      IF(IK.EQ.1) K=2
      IF(IK.EQ.2) K=1
      IMINS=IMINS+1
      GOTO 80001

C      =====
C      CLEAN UP AND EXIT
C      =====

9100   CONTINUE
      TYPE 500,RESET
      CLOSE(UNIT=2)
      CALL QIO("12,7)
      TYPE 91000
91000  FORMAT (////////,10X,'SELECT ONE OF THE FOLLOWING [ENTER 1 OR 2]',
1      /20X,'1 TO SAVE DIGITIZED FILE ',
2      /20X,'2 TO DELETE DIGITIZED FILE ',)
      ACCEPT 617,ANS
      IF(ANS.EQ.'1') THEN
      CLOSE(UNIT=1)
      ELSE IF(ANS.EQ.'2') THEN
      CLOSE(UNIT=1, STATUS='DELETE')
      ELSE
      GO TO 9100
      END IF
      CALL ACSTAT(ISTAT3,JSTAT,1, 5 ,0)
      CALL QIO("12,7)
      CALL QIO("12,9)
      END
      SUBROUTINE DIRECT (TCOUNT,K)
C      STORAGE SEQUENCE AS FOLLOWS:
C      FIRST - RAW RESPIRATION
C      SECOND - DELTA Z
C      THIRD - DZ/DT
C      FOURTH - ECG

```

```

c      TK3 FOR TASK BUILD
C
c      TK3>FILENAME=FILENAME
C      TK3>/
c
      PARAMETER (IFREQ=200)
      integer buf(8*IFREQ), tcount
      COMMON BUF

C      do 100 tcount = 1,2
      WRITE(1,REC=TCOUNT,ERR=999)(BUF((K-1)*4*IFREQ+ICOUNT),
1ICOUNT=1,4*IFREQ)
C          TYPE *, 'OK ON READ'
C          do 200 i = 1, 200
C              write(5,10)i,(BUF(ICOUNT,1),ICOUNT=(i-1)*4+1,i*4)
C200          continue
C          TYPE *, 'OK ON WRITE'
100      continue

10      FORMAT(x,i3,4(2x,i4))
20      FORMAT(1X,I5)

999      CONTINUE
      RETURN
      END
      SUBROUTINE SETFRQ(CFREQ,TFREQ)
C
C*****
C
C      SUBROUTINE GETFRQ(CFREQ,TFREQ)
C
C      LANGUAGE: FORTRAN-77
C
C
C      FUNCTION:
C      This SUBROUTINE will prompt the user for the desired clock frequency to
C      be used in acquiring analog data.
C
C
C      OUTPUTS:
C      CFREQ = REAL*4 variable containing the user's desired frequency in hz.
C      TFREQ = REAL*4 variable containing the user's actual frequency in hz.
C
C
C      SUBROUTINES REFERENCED:
C          SUBROUTINE CLKFRQ.
C
C
C      *****
C      *****
C
C

```

```

REAL*4 CFREQ,TFREQ          !Declare desired, actual frequencies.
C
1    CONTINUE
    CALL CLKFRQ(CFREQ,TFREQ)
    IF(TFREQ.NE.-999.0) GO TO 2      !Skip ahead if CFREQ is ok.
C
    TYPE 9015
9015  FORMAT(/,1X,'Bad frequency, please try again',/)
    GO TO 1          !Prompt for frequency again.
C
2    CONTINUE
999   RETURN          !Return to caller.
    END
    PROGRAM    IMPSTR

C
C    AN IMPEDANCE CARDIOGRAPH PROGRAM USING DIGITIZED DATA FILES
C
    PARAMETER (IFREQ=200)
    DIMENSION SLOPE(2),ENTRCP(2)
C    COMMON ISPACE(4*IFREQ ),SPACE(30*IFREQ )
    REAL*8      DRAW ,IMPACQ ,IMPCAL

C    =====
C    DATA ACQUISITION OF BIOFEEDBACK PARAMETERS
C    =====

    REAL TEMPRT,
@      ZERO,CFREQ,TFREQ,EPOCH

    INTEGER PREDIA,PREHR,PRETEM,PRESYS,PSTDIA,PSTSYS,NUMDIA,
@      PSTHR,PSTTEM,RUNNUM,IDAY,NSN,AGE,SUSC,
@      TESTIM,MIN,PSTTST,PART,CHANEL,TOTMIN,PRETST,
@      IMINS,IAN,          HIGH( 2,2)          ,DIAG(15),
@      CLEANS(4),ICLEAN(2)

    LOGICAL*1 SEX,ANS,TEMP,DMY,IRSP,IABORT,IPAWS,TYMEAN,STDDEV,
@      GROUP(8),DATE(7),DIRECT(2),OSN(2),STRNG(8),
@      FRMT(12)

    CHARACTER*1 CMT(80),RUN(4),RUNTYP,DATFIL(14),EXP
    CHARACTER*14 CATFIL
    EQUIVALENCE (DATFIL,CATFIL)
    DATA CLEANS/27,91,50,74/
C    OPEN (UNIT=1,NAME='NUDRAW.TEL',FORM='FORMATTED',TYPE='OLD')
C    READ(1,2) SLOPE,ENTRCP,C,BSA,CATFIL,MAXREC
C2   FORMAT (6(G15.7),A14,I10)
C    CLOSE(UNIT=1)
9115 CONTINUE
C    TYPE 5020
5020 FORMAT(/' EXPERIMENT DESIGNATION (A,B,...): ', $)
C    ACCEPT 5021,EXP
5021 FORMAT (A1)
    TYPE 5030
5030 FORMAT(/' GROUP TYPE (maximum of 8 characters): ', $)
    ACCEPT 5031,(GROUP(I),I=1,8)

```



```

5031  FORMAT(8A1)
C      TYPE 5040
5040  FORMAT(/' RUN TYPE (maximum of 4 characters): ', $)
C      ACCEPT 5041, (RUN(I), I=1, 4)
5041  FORMAT (4A1)
      TYPE 5050
5050  FORMAT(/' RUN NUMBER: ', $)
      ACCEPT *, RUNNUM
      TYPE 5060
5060  FORMAT(/' DIRECTION (cc,cw,nd): ', $)
      ACCEPT 5061, (DIRECT(I), I=1, 2)
5061  FORMAT (2A1)
C      TYPE 5070
5070  FORMAT(/' DAY NUMBER: ', $)
C      ACCEPT *, IDAY
C      TYPE 5080
5080  FORMAT(/' NEW SUBJECT NUMBER: ', $)
C      ACCEPT *, MSN
      TYPE 5090
5090  FORMAT(/' SUBJECT'S INITIALS (first and last): ', $)
      ACCEPT 5091, (OSN(I), I=1, 2)
5091  FORMAT (2A1)
      TYPE 5100
5100  FORMAT(/' AGE: ', $)
      ACCEPT *, AGE
      TYPE 5110
5110  FORMAT(/' SUSCEPTIBILITY (1,2,3): ', $)
      ACCEPT *, SUSC
      TYPE 5120
5120  FORMAT(/' SEX (m/f): ', $)
      ACCEPT 5121, SEX
5121  FORMAT (A1)
      TYPE 5130
5130  FORMAT(/' TEST DATE (mmnddyy): ', $)
      ACCEPT 5131, (DATE(I), I=1, 7)
5131  FORMAT (7A1)
      TYPE 5140
5140  FORMAT(/' TEST TIME (military in hrs and mins, e.g. for 1:30 p.m.'
      @      ', enter 1330): ', $)
      ACCEPT *, TESTIM
      TYPE 5150
5150  FORMAT(/' PRE HEART RATE (in beats/min): ', $)
      ACCEPT *, PREHR
      TYPE 5160
5160  FORMAT(/' PRE TEMPERATURE: ', $)
      ACCEPT *, TEMPRT
      PRETEM = IIFIX(TEMPRT*10.0)
      TYPE 5170
5170  FORMAT(/' PRE B.P. (sys,dias ,e.g. 120,80): ', $)
      ACCEPT *, PRESYS, PREDIAS
      TYPE 5180
5180  FORMAT(/' PRE/POST TEST BASELINE (in minutes): ', $)
      ACCEPT *, PRETST
      PSTTST = PRETST

```

```

CATFIL(11:14)='.IMP'
CALL STORE (CATFIL,IFLAG,GROUP,RUNNUM,DIRECT,OSN,AGE,SUSC,SEX,
1 DATE,TESTIM,PREHR,TEMPRT,PRETEM,PRESYS,PREDIA,PRETST,PSTTST)
TYPE 500,CLEANS
500 FORMAT(X,4A1)
END
SUBROUTINE STORE (CATFIL,IFLAG,GROUP,RUNNUM,DIRECT,OSN,AGE,
1 SUSC,SEX,DATE,TESTIM,PREHR,TEMPRT,PRETEM,PRESYS,PREDIA,PRETST,
2 PSTTST)

C
C AN IMPEDANCE CARDIOGRAPH PROGRAM USING DIGITIZED DATA FILES
C

DIMENSION SLOPE(2),ENTRCP(2)
COMMON DATCAR(4,2)
REAL*8 DRAW

C
C =====
C DATA ACQUISITION OF BIOFEEDBACK PARAMETERS
C =====

REAL TEMPRT,MEAN(4),SD(4),
@ ZERO,CFREQ,TFREQ,EPOCH

INTEGER PREDIA,PREHR,PRETEM,PRESYS,PSTDIA,PSTSYS,NUMDIA,
@ PSTHR,PSTTEM,RUNNUM,IDAY,NSN,AGE,SUSC,
@ TESTIM,MIN,PSTTST,PART,CHANEL,TOTMIN,PRETST,
@ IMINS,IAN, HIGH( 2,2) ,DIAG(15),
@ CLEANS(4),ICLEAN(2)

LOGICAL*1 SEX,ANS,TEMP,DMY,IRSP,IABORT,IPAWS,TYMEAN,STDDEV,
@ GROUP(8),DATE(7),DIRECT(2),OSN(2),STRNG(8),
@ FRMT(12)
CHARACTER*1 CMT(80),RUNTYP,DATFIL(14),EXP
CHARACTER*4 RUN
CHARACTER*14 CATFIL,IMPFIL
EQUIVALENCE (DATCAR,MEAN),(DATCAR(1,2),SD)

DATA DRAW/'DRAWRM'/ ,ICLEAN/27,99/
DATA CLEANS/27,91,50,74/

C IFLAG=0
C TYPE 500,CLEANS
500 FORMAT(X,4A1)

OPEN (UNIT=1,STATUS='OLD',NAME='DW1:NUDRAW.TEL',FORM='FORMATTED')
READ(1,2) SLOPE,ENTRCP,C,BSA,IMPFIL,MAXREC
2 FORMAT(6( G15.7),A14,I10)
CLOSE(UNIT=1,DISPOSE='SAVE')
IMPFIL(12:14)='IMP'
MIN=MAXREC/4 ! THERE ARE 4 EPOCHS PER MINUTE
N=1 !SET FOR FIRST DISKETTE

```

```

C      -----
C      ASSEMBLE DATAFILE NAME
C      -----

9115  CONTINUE
      CALL ERRSET(63,.TRUE.,.FALSE.,.FALSE.,.FALSE.,MAX)
      CALL ERRSET(36,.TRUE.,.FALSE.,.TRUE.,.FALSE.,MAX)
      CALL ERRSET(29,.TRUE.,.FALSE.,.TRUE.,.FALSE.,MAX)
      EXP=IMPFIL(5:5)
      DECODE(2,9112,IMPFIL(6:7)) NSN
      DECODE(2,9112,IMPFIL(8:9)) IDAY
      IF (IMPFIL(10:10) .EQ. '0') RUN(1:4)='BSLN'
      IF (IMPFIL(10:10) .EQ. '1') RUN(1:4)='CSSI'
      IF (IMPFIL(10:10) .EQ. '2') RUN(1:4)='TRAN'
      IF (IMPFIL(10:10) .EQ. '3') RUN(1:4)='    '
      IF (IMPFIL(10:10) .EQ. '4') RUN(1:4)='TASK'
      IF (IMPFIL(10:10) .EQ. '5') RUN(1:4)='AMBL'
      IF (IMPFIL(10:10) .EQ. '6') RUN(1:4)='    '
      IF (IMPFIL(10:10) .EQ. '7') RUN(1:4)='VARD'
      IF (IMPFIL(10:10) .EQ. '8') RUN(1:4)='DRUM'
      IF (IMPFIL(10:10) .EQ. '9') RUN(1:4)='ZERO'

      CATFIL(1:1) = 'D'
      CATFIL(2:2) = 'Z'
      CATFIL(3:3) = '1'
      CATFIL(4:4) = ':'
      CATFIL(5:10)=IMPFIL(5:10)
9111  FORMAT(A2)
9112  FORMAT( I2)
      CATFIL(11:11) = '.'

C      -----
C      OPEN DATA FILES
C      -----

      TYPE 500,CLEANS
40    TYPE 9113
9113  FORMAT(////////// ' PLEASE INSERT DATA FLOPPY WITH AT
1    LEAST 60 CONTIGUOUS FREE BLOCKS//' INTO DRIVES DZ1 & DZ2.',
1    ' WHEN READY, PRESS RETURN',/////////)
      ACCEPT 10,ITEMP
10    FORMAT(I5)
C
9510  CONTINUE
C
      CATFIL(12:12) = 'I'
      CATFIL(13:13) = 'M'
      CATFIL(14:14) = 'P'
      IREC=0
      CALL ERRSNS
      OPEN (UNIT=1, TYPE='OLD', NAME=CATFIL, FORM='FORMATTED',
1    ERR=45
      )
      GOTO 100
45    CONTINUE

```

CHAPTER 10
OF POOR QUALITY

```

CALL ERRSNS (IERNUM)
IF(IERNUM.EQ.29) THEN
  OPEN (UNIT=1, TYPE='NEW', NAME=CATFIL, FORM='FORMATTED',
1 ERR=40 )
  ELSE
  GOTO 40
END IF
100 CONTINUE
C
C -----
C STORE KEY FIELD INFO TO DATA FILE
C -----
D TYPE *, 'EPOCHS, RUN DURATION (MINUTES)= ',IMINS,MIN
WRITE(1,9120) EXP,(GROUP(I),I=1,8),(RUN
1 RUNNUM,(DIRECT(I),I=1,2),IDAY,NSN,(OSN(I),I=1,2),
2 AGE,SUSC,SEX,(DATE(I),I=1,7),TESTIM,PREHR,PRETEM,PRESYS,PREDIA,
3 PSTHR,PSTTEM,PSTSYS,PSTDIA,MIN,PRETST,
4 MIN-PRETST-PSTTST,PSTTST,NUMDIA, 1 , 8 ,DIAG(1)

9120 FORMAT( A1,X,8A1,X, A4,X,I2,X,2A1,2(X,I2),X,2A1,X,I2,X,I1,
1 X,A1,X,7A1,X,I4,X,I3,X,I4,X,I3,X,I3,
2 X,I3,X,I4,5(X,I3).5(X,I2))

DO 9130 I = 2,14
  WRITE(1,9140) DIAG(I)
9130 CONTINUE
9140 FORMAT(109X,I2)

WRITE(1,9150) DIAG(15),MAXREC
9150 FORMAT(109X,I2, I3)

FRMT(2) = '0'
FRMT(3) = '6'
ZERO = 0.0
IERR=0
9650 CONTINUE
IREC=IREC+1
IF(IREC.GT.MAXREC) GOTO 9165
D TYPE *, ' IREC= ',IREC,' MAXREC= ',MAXREC
CALL DRECT (IREC,IERR,IMPFIL)
IF(IERR.EQ.1) GOTO 9160
WRITE(1,9550) MEAN,0.,0.,0.,0.,SD,0.,0.,0.,0.
C WRITE(1,9550) (DATCAR(I,1),I=1,4),0.,0.,0.,0.,
C 1 (DATCAR(I,2),I=1,4),0.,0.,0.,0.
9550 FORMAT(8(X,F9.4))
GOTO 9650
9160 CONTINUE
WRITE(1,9550) 0.,0.,0.,0.,0.,0.,0.,0.,0.,0.,0.,0.,0.,0.,0.
GOTO 9650
9165 CONTINUE
CLOSE(UNIT=1)
CATFIL(12:12)='C'
C DATFIL(13)='A'

```

```

C      DATFIL(14)='L'
C      OPEN (UNIT=1, TYPE='NEW', NAME=CATFIL, FORM='FORMATTED',
C      1 ERR=40
C      DO 9700 I=1,3
C      WRITE(1,9600) HIGH(1.2),HIGH(1.1),CALVAL(1,I),CALVAL(2,I)
C9600   FORMAT( 2(X, I10),2(X,F10.4))
C9700   CONTINUE
C      CLOSE(UNIT=1)
C      IF(IFLAG.EQ.0) GOTO 9875
C      CATFIL(13:13)='O'
C      CATFIL(14:14)='M'
C      OPEN (UNIT=1, TYPE='NEW', NAME=CATFIL, FORM='FORMATTED',
C      1 ERR=40
C      OPEN (UNIT=2, TYPE='OLD', NAME='SY:TEMP.COM', FORM='FORMATTED',
C      1 READONLY
C      REWIND 2
C      J=0
C      KCNT=80
9750   CONTINUE
C      J=J-KCNT
C      IF(J.GT.KCHAR) KCNT=KCHAR-J+KCNT
C      READ(2,9800) (CMT(I),I=1,KCNT)
C      WRITE(1,9801)(CMT(I),I=1,KCNT)
9800   FORMAT(80A1)
9801   FORMAT(X,80A1)
C      IF(J.GE.KCHAR) GOTO 9850
C      GOTO 9750
9850   CONTINUE
C      CLOSE(UNIT=1)
C      CLOSE(UNIT=2)
9875   CONTINUE
C      IF(N.EQ.2) GOTO 99990
C      N=2
C      CATFIL(3:3)='2'
C      GOTO 9510
99990   CONTINUE
C      TYPE *, ' ARE YOU FINISHED WITH THE DIGITIZED DATA ?[Y/N] '
C      ACCEPT 10000,ANS
10000   FORMAT(A1)
C      IF(ANS.EQ.'Y') THEN
C      TYPE *, ' THEN YOU WANT TO DELETE THE DIGITIZED FILE ?[Y/N]'
C      ACCEPT 10000,ANS
C      IF(ANS.EQ.'Y') THEN
C      OPEN (UNIT=2,NAME=IMPFIL, TYPE='OLD')
C      CLOSE (UNIT=2,DISPOSE='DELETE')
C      IMPFIL(11:14)='.DIG'
C      OPEN (UNIT=2,NAME=IMPFIL, TYPE='OLD')
C      CLOSE (UNIT=2,DISPOSE='DELETE')
C      END IF
C      END IF
C      IF(IFLAG.EQ.1) THEN
C      OPEN (UNIT=2,NAME='SY:TEMP.COM',TYPE='OLD')
C      CLOSE (UNIT=2,DISPOSE='DELETE')
C      END IF

```

```

C      OPEN (UNIT=1,STATUS='OLD',NAME='DWI:NUDRAW.TEL',FORM='FORMATTED')
C      CLOSE(UNIT=1,DISPOSE='DELETE')
C      TYPE 91100
91100  FORMAT (/' SUCCESSFUL COMPLETION ')
99999  CONTINUE
      RETURN
      END
      SUBROUTINE DRECT (TCOUNT,IERR,IMPFIL)
C      STORAGE SEQUENCE AS FOLLOWS:
C      FIRST - RAW RESPIRATION
C      SECOND - DELTA Z
C      THIRD - DE/DT
C      FOURTH - ECG

c      TKB FOR TASK BUILD
C
c      TKB>FILENAME=FILENAME
C      TKB>/
C      TKB>MAXBUF=1600
C      TKB>//
c

      PARAMETER (IFREQ=200)
C      integer buf(4*IFREQ,1)
      INTEGER tcount
C      DIMENSION X(15*IFREQ),Y(15*IFREQ)
      DIMENSION DATCAR(4,2)
      COMMON DATCAR
C      EQUIVALENCE (DATCAR,BUF)
      CHARACTER*14 IMPFIL

      OPEN (UNIT=2, status='old', NAME=IMPFIL
1      FORM='UNFORMATTED', access='direct',recl=8 )

C      do 100 tcount = 1,2
          read(2,REC=TCOUNT,ERR=990)((DATCAR(I,J),I=1,4),J=1,2)
C          TYPE *, 'OK ON READ'
C          do 200 i = 1, 200
C              write(5,10)i,(BUF(ICOUNT,1),ICOUNT=(i-1)*4+1,i*4)
C200      continue
C          TYPE *, 'OK ON WRITE'
100      continue
      GOTO 999
990      CONTINUE
      IERR=1
999      CLOSE(UNIT=2, STATUS='SAVE')
10      FORMAT(x,i3,4(2x,i4))
20      FORMAT(1X,I5)

      RETURN
      END
      PROGRAM    DRAWED
      PARAMETER (IFREQ=200) ! SET SAMPLING SIZE FOR DIGITIZED DATA
      DIMENSION ISTAT(2),IRES(75),OUTPUT(10,70),

```

```

1  SLOPE(2),ENTRCP(2)
   CHARACTER*1 ANS
   CHARACTER*1 ICHRPK(2),IBUF,ISTOP,IREFET
   CHARACTER*4 EPOCH,ETIME
   VIRTUAL IDATA(9000)
   DIMENSION Y(1450),X(1000),NCHAN(3),SMIN(3),SMAX(3),DATCAR(5,50),
1  DATSD(4),IPAR(6)
   INTEGER*2 IRESP(IFREQ)
   CHARACTER*14 CATFIL
   INTEGER LPTS(3),IPEAKS(50),IOSB(2),CLEANS(2)
   BYTE I233,I73,I162,IHOME(4)
   CHARACTER*2 IPT,IPB
   COMMON /PLOT/X,Y
   EQUIVALENCE (DATCAR,Y(1201)),(Y(501),OUTPUT)
   DATA NCHAN/3,4,1/,SMIN/350.,150.,0./,SMAX/550.,350.,150./
   DATA LPTS/3,3,10/,NSMFRQ/100/,NREC/15/,IREAD/"001000/
   DATA A/.01/,B/0./,C/0./,IBRAC/' '],ICHAN/3/
   DATA IGATE/0/,IM/0/,I233/"233/,I73/"73/,I162/"162/
   DATA IHOME/155,63,54,108/,CLEANS/27.99/
C  CALL GETADR(IPAR(1),IBUF) !GET START ADDR OF IBUF & STOR IN IPAR(1)
C  IPAR(2)=1 !SET TO BYTE SIZE OF IBUF
   TYPE 4000,CLEANS
4000  FORMAT(X,4A1)
   TYPE 4500
4500  FORMAT (10X,60('*'),/////,22X,' DATA REDUCTION & EDITING ',
1  ////,10X,60('*'),////)
C
C  READ FILE CONTAINING SLOPES, INTERCEPTS, C AND BSA
C      WHERE  C =  $\rho \cdot L^{**2}$ 
C              BSA = BODY SURFACE AREA
C
   OPEN (UNIT=1,NAME='DW1:NUDRAW.TEL',FORM='FORMATTED',TYPE='OLD')
   READ(1,2) SLOPE,ENTRCP,C,BSA,CATFIL,MAXREC
2  FORMAT(6( G15.7),A14,I10)
   CLOSE (UNIT=1)

   IDZDT0=-ENTRCP(2)/SLOPE(2)! GET AD VALUE FOR DZDT=0
   IDZDT1=(1.-ENTRCP(2))/SLOPE(2)! GET AD VALUE FOR DZDT=1
100  CONTINUE

   TYPE *, ' STARTING TIME CODE FOR DATA ON TAPE (eg. enter 300 for',
1  ' 3:00)'
   ACCEPT *,ISTREC
   TYPE *, ' EPOCH FOR START OF DATA REDUCTION '
   ACCEPT *,ISTEPC
   TYPE *, ' ENTER LAST DATA EPOCH '
   ACCEPT *,IENEPC

C
C  KREC SETS THE STARTING RECORD FOR DATA IN DIGITIZED FILE
C  INDREC SETS THE ENDING RECORD FOR DATA IN DIGITIZED FILE
C  IDKREC SETS THE 15 SECOND EPOCH NUMBER (DISPLAYED TO SCREEN)
C      AND USED TO STORE RESULTS IN 15 SEC EPOCHS
C

```

```

KREC= (ISTEPC-1)*15+61
INDREC=(IENEPC)*15+61
MAXREC=(INDREC-61)/15
IDKREC=ISTEPC
CATFIL(11:14)='.DIG'

105  CONTINUE
C    CALL QIO(IREAD,7,4,,IOSB,IPAR,IDSW)
10   CONTINUE
      XMIN=0.
      XMAX=IFREQ *NREC
      CALL CGL( 90) !INITIALIZE GRAPHICS
      CALL CGL( 92) !NEW FRAME
      CALL CGL(103,'DW1:FILE2.GID',13) !INITIALIZE PLOT FILE
      CALL CGL(105,'DW1:FILE2.GID',13) !SELECT PLOT FILE
      CALL CGL( 80,XMIN,XMAX,0.      ,600.  ) !SET WINDOW TO DEFAULT VALUES
      CALL CGL(86,0) !SET ORIGIN
      OPEN (UNIT=' ', status='old', NAME=CATFIL, !OPEN DIGITIZED DATA FILE
1      FORM='UNFORMATTED', access='direct',recl=4*IFREQ)
C    DO 125 I=1,NREC
C      IREC=KREC+I-1
C
C      read(1,REC=IREC,ERR=999)((IDATA(J+(I-1)*IFREQ),ID,ID,ID),
C      1  J=1,IFREQ)
125  CONTINUE
C998  CLOSE(UNIT=1, STATUS='SAVE')
C      M=0
C      NPTS=LPTS(3)
C      DO 135 L=1,IFREQ*NREC,NPTS
C        M=M+1
C        X(M)=L
C        Y(M)=IDATA(L)
C135  CONTINUE
      M=0 !SET COUNTER OF RESPIRATION VALUES
      IFRQRC=IFREQ*NREC !SET AREA SIZE FOR DIG (IDATA)BLOCKS
      DO 150 I=1,NREC
      IREC=KREC+I-1

C      STORAGE SEQUENCE OF DIGITIZED DATA AS FOLLOWS:
C      FIRST - RAW RESPIRATION - STORED IN RESP FROM 1 TO IFREQ
C      SECOND - DELTA Z - STORED IN IDATA FROM 1 TO IFRQRC
C      THIRD - DZ/DT - STORED IN IDATA FROM 1+IFRQRC TO 2*IFRQRC
C      FOURTH - ECG - STORED IN IDATA FROM 1+2*IFRQRC TO 3*IFRQRC

      read(1,REC=IREC ,ERR=999)((IRESP(J)
1 IDATA(J+(I-1)*IFREQ),IDATA(IFRQRC+J+(I-1)*IFREQ),
2 IDATA(2*IFRQRC+ J+(I-1)*IFREQ)),J=1,IFREQ)

C      PRINT *,(IECG(J+(I-1)*IFREQ ),J=1,IFREQ),J,I
C
C      STORE A SUBSET( EVERY LPTS(3) VALUE) OF RESP IN Y FOR PLOTTING
C
      DO 175 J=1,IFREQ,LPTS(3)
      M=M+1

```



```

      Y(M)=IRESP(J)
175  CONTINUE
150  CONTINUE
999  CLOSE(UNIT=1, STATUS='SAVE') !CLOSE DIG. DATA FILE
C
C   PLOT FIRST RESPIRATION(K=3) THEN ECG(K=2) AND LAST DZ/DT(K=1)
C
      DO 3000 K=3,1,-1
      YSMIN=SMIN(K) !PARTITION SCREEN AREA MIN FOR K PLOT
      YSMAX=SMAX(K) !PARTITION SCREEN AREA MAX FOR K PLOT
      NPTS=LPTS(K)
C
C   -----
C   SET PLOT ARRAY PAIRS (X,Y)
C   -----
      IF(K.EQ.3) THEN
      DO 275 L=1,M !SETUP X VALUES FOR RESPIRATION PLOTTING
      X(L)=1+(L-1)*NPTS
275  CONTINUE
      ELSE
      M=0
      DO 300 L=1,IFREQ *NREC,NPTS !SETUP (X,Y) PAIRS FOR ECG AND DZ/DT PLOT
      M=M+1
      X(M)=L
      Y(M)=IDATA(IFRQRC*K+L)
300  CONTINUE
      END IF
      CALL INDEX(XMIN,XMAX,YMIN,YMAX,M) !DETERMINE PLOT LIMITS
      IF(K.EQ.1) THEN
      IF(IDZDT1-IDZDTO.LT.(YMAX-YMIN)/3) THEN !REDONE LIMITS IF NEEDED
      YMAX=IDZDT1+IDZDT1-IDZDTO
      YMIN=IDZDTO-(IDZDT1-IDZDTO)
      END IF
      IF(IDZDT1.GT.YMAX) YMAX=IDZDT1
      DZDTO=(IDZDTO-YMIN)/(YMAX-YMIN)*(YSMAX-YSMIN)+YSMIN !SCALE DZDTO
      DZDT1=(IDZDT1-YMIN)/(YMAX-YMIN)*(YSMAX-YSMIN)+YSMIN !SCALE DZDT1
C
C   PLOT DZ/DT CALIBRATION VALUES FOR OBSERVATION
C
      CALL CGL( 1,X(1), DZDTO) !MOVE "PEN" TO POSITION
      CALL CGL( 4,XMAX, DZDTO) !DRAW TO POSITION
      CALL CGL( 1,X(1), DZDT1) !MOVE "PEN" TO POSITION
      CALL CGL( 4,XMAX, DZDT1) !DRAW TO POSITION
      END IF
C
      DO 322 L=1,10
C
      TYPE *,X(L),Y(L),YMIN,YMAX,L
C
      Y(L)=(Y(L)-YMIN)/(YMAX-YMIN)*(YSMAX-YSMIN)+YSMIN
C
      TYPE *,X(L),Y(L),YMIN,YMAX,L
C322  CONTINUE
      DO 325 L=1,M !SCALE Y
C
      PRINT *,X(L),Y(L),YMIN,YMAX,L
      Y(L)=(Y(L)-YMIN)/(YMAX-YMIN)*(YSMAX-YSMIN)+YSMIN
C
      TYPE *,X(L),Y(L),YMIN,YMAX,L
325  CONTINUE
C
      FACT=5./32767.

```

```

      CALL CGL(1,X(1),Y(1)) MOVE PEN TO POSITION
C      KOLD=X(600)+1./200.
C      NY=NREC*IFREQ
      NYPLT=IFRQRC /NPTS
      CALL CGL(6,X,Y,NYPLT) !PLOT NYPLT X,Y PAIRS
      IF (K.EQ.1) THEN !PLOT EPOCH NUMBER AND START TIME OF EPOCH
      IEPOCH=IMOD(KREC-61,60)+(KREC-61)/60*100+ISTREC
      ENCODE (4,32500,ETIME)IEPOCH
      ENCODE (4,32500,EPOCH)IDKREC
32500  FORMAT(I4)
      CALL CGL (1,XMIN,SMIN(1))
      CALL CGL (16,EPOCH,4)
      CALL CGL (1,XMIN+240.,SMIN(1))
      CALL CGL (16,ETIME,4)
      END IF
500    CONTINUE
3000   CONTINUE !END OF PLOTTING
C
C      NEXT CALL PEAK DETECTION
C
D      OPEN (UNIT=2,STATUS='NEW',NAME='DW1:IMPTEMP.TST',
D      1  FORM='FORMATTED',ACCESS='SEQUENTIAL')
      CALL VOL(IFRQRC,NPEAKS,IDATA,IADMIN,IGATE,IM,IDZDTO,IDZDT1)
C      IDZDT - INPUT OF DZ/DT DATA
C      IFRQRC - NUMBER OF DZ/DT DATA POINTS
C      OUTPUT(3,I) - ARRAY CONTAINING TIME VALUES FOR DETECTED PEAKS
C      NPEAKS - NUMBER OF PEAKS DETECTED
      IF(NPEAKS.GT.0)THEN
      CALL CGL(12,4,0,0) !SET LINE STYLE TO DOTTED
C
C      PLOT PEAK INDICATORS
C
      DO 680 I=1,NPEAKS
      YLOW=SMIN(3)+40
      XL=(OUTPUT(3,I)-1.)
      CALL CGL(1,XL,SMAX(1) ) !POSITION TO X=XL,Y=SMAX(1)
      CALL CGL(4, XL,YLOW) !DRAW LINE TO X=XL,Y=YLOW
      ENCODE (2,610,ICHRPK),I
610    FORMAT( I2)
      XLOW=XL-10.
      DO 620 J=1,2
      CALL CGL(1,XLOW,YLOW )
      CALL CGL(16,ICHRPK(J),1) !WRITE PEAK NUMBER J
      YLOW=YLOW-20.
620    CONTINUE
680    CONTINUE
      CALL CGL(12,1,0,0) !RESET LINE STLYE TO SOLID
C
C      USE A SUBROUTINE HERE TO CALCULATE THE TIME WHEN
C      dZ/dT=0 BEFORE EACH PEAK ALSO CALCULATE STROKE VOLUME
C
      CALL TFACT(NPEAKS,IFREQ ,C,NY      ,K,TYMEAN,
1      SLOPE,ENTRCP,BSA,IDATA,KREC,IDZDTO,IADMIN,CATFIL,IDKREC)
C      ELSE

```

```

C      SET MEAN AND SD TO ZERO
      END IF
C      IF(ICOUNT.GT.1) THEN
C      CALL CGL(12,4,0,0)
C      DO 700 I=1,ICOUNT-1
C      XL=600.*I/200.
C      CALL CGL(1,XL,YMAX      )
C      CALL CGL(4,XL      ,YMIN      )
C700   CONTINUE
D      CLOSE(UNIT=2)
C
C      SET UP ONE LINE DIALOG AREA NOW
C
      L20=49
      L24=49
      WRITE(5,40000) I233,L20,I73,L24,I162
40000  FORMAT(1X,5A1)
      DO 725 I=1,50 !ZERO ALL IPEAKS - ARRAY OF USER SELECTED PEAKS
      IPEAKS(I)=0
725   CONTINUE
C
C      QUIRY EVENT FLAG 4 TO SEE IF USER REQUEST MADE
C
C      CALL READEP(4,IEF)
C      IF(IEF.GT.0) THEN ! REQUEST MADE VIA KEYBOARD
C      END IF
C      CALL QIO("12,7) ! CANCEL QUE I/O REQUEST
C
C      DIALOG SECTION
C
      TYPE *, ' SKIP THIS EPOCH ? [ENTER Y/N]'
      ACCEPT 10000,ANS
      IF(ANS.EQ.'Y') THEN
      KREC=KREC+NREC !INCREMENT RECORD INDEX
      MAXREC=MAXREC-1 !REDUCE MAX RECORD SIZE BY ONE
      CALL CGL(106,'DW1:FILE2.GID',13) !DESELECT PLOT FILE
      CALL CGL(104,'DW1:FILE2.GID',13) !TERMINATE PLOT FILE
      CLOSE (UNIT=1,DISPOSE='DELETE')
      IF (KREC.GE.INDREC) GOTO 99999
      GOTO 105
      END IF
      TYPE *, ' REDO ANALYSIS ? [ENTER Y/N]'
      ACCEPT 10000,ANS
      IF(ANS.EQ.'Y') THEN
      TYPE *, ' ENTER GATE FACTOR, MINIMUM DIFFERENCE FACTOR NOW'
      ACCEPT *,IGATE,IM
      CALL CGL(106,'DW1:FILE2.GID',13)
      CALL CGL(104,'DW1:FILE2.GID',13)
      CLOSE (UNIT=1,DISPOSE='DELETE')
      GOTO 10
      END IF
      TYPE *, ' OK TO CALCULATE HEATHER INDEX ? [Y/N]'
      ACCEPT 10000,ANS

```

ORIGINAL PAGE IS
OF POOR QUALITY

```

IF(ANS.EQ.'Y') THEN
CALL HICALC(NPEAKS,IFREQ ,IDATA ,KREC,NREC,CATFIL)
ELSE
DO 7250 I=1,NPEAKS
7250 DATCAR(4,I)=0.0
CONTINUE
END IF
CALL CGL(106,'DW1:FILE2.GID',13)
CALL CGL(104,'DW1:FILE2.GID',13)
CALL CGL(93,INAM,ICODE) !REPORT ERROR
IF(ICODE.NE.0) THEN
TYPE *, ' INAM= ',INAM,' ICODE= ',ICODE
CALL CGL( 91) !TERMINATE GRAPHICS
END IF
7025 CONTINUE
CALL ERRSNS
C PRINT *, ' NPEAKS= ',NPEAKS
C PRINT *, (IPEAKS(II),II=1,NPEAKS)
TYPE *, ' LIST SELECTED PEAKS NOW (SEPERATE W/ COMMA )'
READ(5,7725) (IPEAKS(I),I=1,NPEAKS)
7725 FORMAT(<NPEAKS>I3)
730 CONTINUE
740 CONTINUE
C PRINT *, ' NPEAKS= ',NPEAKS
C PRINT *, (IPEAKS(II),II=1,NPEAKS)
KPEAKS=0
DO 750 I=1,NPEAKS !DETERMINE NUMBER OF PEAKS SELECTED
IF(IPEAKS(I).EQ.0) GO TO 775
KPEAKS=KPEAKS+1
750 CONTINUE
775 CONTINUE
C TYPE *, ' KPEAKS= ',KPEAKS ,(IPEAKS(I),I=1,KPEAKS)
IF(KPEAKS.GT.0) THEN
C PRINT *,IDKREC,( IPEAKS(I),(DATCAR(K,IPEAKS(I)),K=1,4),I=1,KPEAKS)
CALL STRVOL(KPEAKS,DATSD,IPEAKS)
ELSE !ZERO STROKE VALUE,CARDIAC OUTPUT,CARDIAC INDEX, HEATHER INDEX
DO 810 I=1,4
DATCAR(I,1)=0.0 !MEANS
DATSD(I)=0.0 !STANDARD DEVIATIONS
810 CONTINUE
END IF
CALL ERRSNS(IERR)
IF(IERR.NE.0) GOTO 7025
820 CONTINUE
C
C WRITE RESULTS OF EPOCH(=IDKREC) CALCS TO INTERNAL FILE
C
CATFIL(11:14)='.IMP'
OPEN (UNIT=1,STATUS='OLD',NAME=CATFIL,
1 FORM='UNFORMATTED',ACCESS='DIRECT',RECL=8 )
WRITE(1,REC=IDKREC) ((DATCAR(K,1),K=1,4)),DATSD
CATFIL(11:14)='.DIG'
CLOSE(UNIT=1)
CLOSE(UNIT=2)

```

```

TYPE *, ' DO YOU WANT A COPY ? [Y/N]'
ACCEPT 10000,ANS
10000 FORMAT(A1)
      IF (ANS.EQ.'Y') THEN
        IRES(1)=1
        CALL CPRNT(ISTAT,IRES,'DW1:FILE2.GID',13)
      ELSE
        OPEN (UNIT=1,STATUS='OLD',NAME='DW1:FILE2.GID')
        CLOSE (UNIT=1,DISPOSE='DELETE')
      END IF
      TYPE *, ' EXIT ? [Y/N]'
      ACCEPT 10000,ANS
      IF(ANS.EQ.'Y') GOTO 99999
      IDKREC=IDKREC+1 !INDEX EPOCH COUNTER FOR PLOT AND RESULTS FILE
      KREC=KREC+NREC !INDEX DIGITIZED DATA FILE RECORD COUNTER
      IF(KREC.GE.INDREC) GOTO 99999 !AT END OF DATA FILE YET
      GOTO 105
99999 CONTINUE
      TYPE 4000,CLEANS
      OPEN (UNIT=1,NAME='DW1:NUDRAW.TEL',FORM='FORMATTED',TYPE='OLD')
      WRITE(1,2) SLOPE,ENTRCP,C,BSA,CATFIL,MAXREC !RECORD MAXREC VALUE
      CLOSE (UNIT=1)
C      WRITE(5,40000) I233,I73,I162
C      TYPE *, ' EXST CALLED ', ' ISTAT= ', ISTAT
      CALL EXST(ISTAT) !SET FLAG FOR EXIT TO PARENT TASK
      STOP
      END
      SUBROUTINE INDEX (XMIN,XMAX,YMIN,YMAX,NPTS)
      PARAMETER (IFREQ=200)
      DIMENSION X(1000),Y(1450)
      COMMON /PLOT/X,Y
C      PRINT *,K,NPTS,(Y(I),I=1,NPTS)
      SX=639./NPTS
      YMIN=Y(1)
      YMAX=YMIN
C      DO 300 I=1,NPTS
      PRINT *,Y(I),I
      IF(Y(I).LT.YMIN) YMIN=Y(I)
      IF(Y(I).GT.YMAX) YMAX=Y(I)
300 CONTINUE
C      SY=479./(IYMAX-IYMIN)
C      ENY=-479.*IYMIN/(IYMAX-IYMIN)
      XMIN=0.
      XMAX=X(NPTS)
C      PRINT *,Y,YMIN,YMAX,NPTS
      RETURN
      END
      SUBROUTINE VOL(NY,KPEAKS,IY,IADMIN,IGATE,IM,IDZDTO,IDZDT1)
C      DO PEAK DETECTION USING DIGITAL ROUTINE 'PEAK'
      PARAMETER (IFREQ=200)
      VIRTUAL IY (9000)
      DIMENSION INPUT(15*IFREQ)
      DIMENSION OUTPUT(10,70)
      DIMENSION ITABLE(68),VTYPE(2,2)

```

ORIGINAL PAGE IS
OF POOR QUALITY

```

      INTEGER BUF( IFREQ,1)
      DIMENSION Y(1450)
      COMMON /PLOT/IX(2000), Y
      EQUIVALENCE (IX,INPUT),(OUTPUT,Y(501))
      DATA VTYPE/' VA','LLEY','BASE','LINE'/

      DATA ITABLE/1, 5, 5,800 ,1,0,0,61*0/
      DATA INLAST,IDIMO/600,70/
      CALL ERRSET(112,.TRUE.,.FALSE.,.FALSE.,.FALSE.,)
      DZREF=IDZDT0+(IDZDT1-IDZDT0)*.5
C      TYPE *, ' HERE IN VOL ',NY,IY
C
C      SET TABLE VALUES FOR DIGITAL PEAK DETECTION ROUTINE
C
C      ITABLE(6)=0
C      ITABLE(7)=0
C      INPTR=0
C      INLAST=NY
C      NPEAKS=0
C      ITABLE(3)=2 !SET GATE FACTOR
C      ITABLE(4)=100 !SET MIN DIFFERENCE FACTOR
C      IF(IGATE.NE.0) THEN !RESET IGATE AND IM IF REQUESTED BY USER
C      ITABLE(3)=IGATE
C      ITABLE(4)=IM
C      END IF
C
C      TRANSCRIBE DZDT DATA TO ARRAY INPUT
C
C      N=0
C      DO 1000 I=1,NY
C      INPUT(I)=IY(3000+I)
C      TYPE *,I,INPUT(I),IY(I)
C 1000    CONTINUE
C      IADMIN=-16000
C      IADMIN=IMIN(INPUT,NY )
C      OPEN (UNIT=2,TYPE='NEW',NAME='INPUT.DAT',FORM='FORMATTED')
C      WRITE (2,1980) (INPUT(I),I=1,NY)
C 1980    FORMAT(8I10)
C      CLOSE (UNIT=2)
C      IF(IADMIN.LT.0) THEN !OFFSET INPUT VALUES SO ALL ARE NON NEGATIVE
C      DO 2000 I=1,NY
C      INPUT(I)=INPUT(I)-IADMIN
C 2000    CONTINUE
C      END IF
C      TYPE *,      (INPUT (I )      ,I=1,NY      ) ,IADMIN
C      WRITE(2,*) (( INPUT (I ) ) ,I=1,NY ) ,IADMIN,SLOPE,INTRCP
C
C      CALL ROUTINE TO FIND PEAKS AND TRAILING MINS TIMES
C      FOR dZ/dT MIN=PEAKS HEIGHT AND T=TIME AT TRAILING
C      MINS - TIME WHEN dZ/dT=0 BEFORE PEAKS
C      OUTPUT ARRAY CONTAINS RESULTS OF PEAK DETECTION
C
C      CALL PEAK(ITABLE,INPUT,INLAST,INPTR,OUTPUT,IDIMO,NPEAKS)
C      TYPE *,NPEAKS,INLAST,INPTR

```

ORIGINAL PAGE IS
OF POOR QUALITY

```

C
C
C      NEXT SECTION NEEDED BECAUSE PEAK ROUTINE MAY AVERAGE INPUT FOR PEAKS

      DO 2500 I=1,NPEAKS
      OUTPUT(2,I)=INPUT(IFIX(OUTPUT(3,I)))
2500  CONTINUE

      IF(IADMIN.LT.0) THEN !CONVERT RESULTS IF DATA OFFSET
      DO 3000 I=1,NPEAKS
      OUTPUT(2,I)=OUTPUT(2,I)+IADMIN
      OUTPUT(4,I)=OUTPUT(4,I)+IADMIN
      OUTPUT(7,I)=OUTPUT(7,I)+IADMIN
3000  CONTINUE
      END IF

C
C      REFINE PEAKS BASED ON SOME CRITERION IF NECESSARY
C

      KPEAKS=0
      DO 4000 I=1,NPEAKS
      IF(OUTPUT(2,I).LE.DZREF) GOTO 4000
      KPEAKS=KPEAKS +1
      IF(KPEAKS.EQ.I) GOTO 4000
      DO 3500 J=1,10
      OUTPUT(J,KPEAKS)=OUTPUT(J,I)
3500  CONTINUE
4000  CONTINUE
D      WRITE (2, * ) NPEAKS,INLAST,INPTR
D      WRITE (2,20000)
C      WRITE (1 ,20000)
20000  FORMAT(' PEAK NO.',8X,'AREA',4X,'P HEIGHT',6X,'P TIME',4X,
A      'L HEIGHT',6X,'L TIME',/,11X,'HALF WIDTH',4X,'T HEIGHT',6X,
B      'T TIME',8X,'TYPE',8X,'RATE'//)
      DO 4 L=1,NPEAKS
      KK=OUTPUT(9,L)+1
D      WRITE(2, 30000)(L,(OUTPUT(I,L),I=1,8),(VTYPE(M,KK),M=1,2),
D      1 OUTPUT(10,L))
C      WRITE(1 , 30000)(L,(OUTPUT(I,L),I=1,8),(VTYPE(M,KK),M=1,2),
C      1 OUTPUT(10,L))
30000  FORMAT(19,5F12.0,/,9X,3F12.0,4X,2A4,F12.0)
4      CONTINUE
      RETURN
      END
      FUNCTION IMIN (INPUT,ICOUNT)
      DIMENSION INPUT(ICOUNT)
      IMIN=INPUT(1)
      DO 1000 I=2,ICOUNT
      IF(INPUT(I).LT. IMIN) THEN
      IMIN=INPUT(I)
      END IF
1000  CONTINUE
      RETURN
      END
      SUBROUTINE TFACT (NPEAKS,IFREQ,C,INLAST,K,TYMEAN,
1  SLOPE,INTRCP,BSA,IZO,IREC,IDZDTO,IADMIN,CATFIL,IDKREC)

```

```

PARAMETER (JFREQ=200)
PARAMETER (NREC=15 )
PARAMETER (NUM=220)
DIMENSION IEPC(NUM),ZOEPC(NUM)
DIMENSION INPUT(15*JFREQ),OUTPUT(10,70),DATCAR(5,50)
VIRTUAL IZO(9000 )
COMMON /PLOT/X(1000),Y(1450)
REAL SLOPE(2),INTRCP(2)
CHARACTER*14 CATFIL
EQUIVALENCE(Y(501),OUTPUT),(X,INPUT),(DATCAR,Y(1201))
LOGICAL*1 TYMEAN
DATA ICOUNT/0/
DATA IEPC/ 1,2,3,4,
15,6,7,8,9,10,11,12,13,14,15,16,17,18,19,20,21,22,23,24,
125,26,27,28,29,30,31,32,33,34,35,36,37,38,39,40,41,42,43,44,45,
146,47,48,49,50,51,52,53,54,55,56,57,58,59,60,61,62,63,64,65,66,
167,68,69,70,71,72,73,74,75,76,77,78,79,80,81,82,83,84,85,86,
1 87, 88, 89, 90, 91, 92, 93, 94, 95 ,96, 97, 98, 99,100,101,102,
1103,104,105,106,107,108,109,110,111,112,113,114,115,116,117,118,
1119,120,121,122,123,124,125,126,127,128,129,130,131,132,133,134,
2135,136,137,138,139,140,141,142,143,144,145,146,147,148,149,150,
3151,152,153,154,155,156,157,158,159,160,161,162,163,164,165,
4166,167,168,169,170,171,172,173,174,175,176,177,178,179,180,
5181,182,183,184,185,186,187,188,189,190,191,192,193,194,195,
6196,197,198,199,200,201,202,203,204,205,206,207,208,209,210,
7211,212,213,214,215,216,217,218,219,220/
DATA ZOEPC/20.9,21.0,21.0,21.0,20.9,21.0,21.0,20.9,20.9,20.9,21.0,
1 21.0,21.0,21.0,21.0,21.0,21.0,21.0,21.0,21.0,21.0,20.9,21.0,
1 21.0,20.9,20.9,21.0,21.0,20.9,20.9,20.9,20.9,20.9,20.9,20.9,20.9,
2 20.9,21.0,21.0,21.1,22.4,22.4,21.3,21.1,21.2,21.2,21.2,21.2,21.2,
3 21.3,21.2,21.2,21.2,21.2,21.5,21.3,21.2,21.3,21.4,21.3,21.2,21.3,
4 21.3,21.3,21.3,21.4,21.3,21.5,22.2,22.2,21.4,21.1,21.3,21.2,21.1,
5 21.1,21.2,21.2,21.2,21.2,21.3,21.2,21.3,21.3,21.3,21.2,21.3,21.3,
6 21.3,21.4,21.4,21.3,21.3,21.3,21.4,21.2,21.2,21.1,21.1,21.1,21.0,
1 21.0,
7 21.1,21.0,21.2,21.3,21.3,21.3,21.3,21.4,21.2,21.3,21.5,21.3,21.4,
8 21.3,21.4,21.5,21.4,21.4,21.5,21.5,21.3,21.2,21.2,21.1,21.2,21.1,
9 21.1,21.1,21.1,21.1,21.4,21.4,21.4,21.3,21.5,21.3,21.3,21.3,21.5,
1 21.4,
1 21.4,21.4,21.5,21.3,21.3,21.3,21.3,21.7,21.3,21.1,21.1,21.7,21.9,
1 21.3,21.3,21.1,21.8,21.1,21.2,21.3,21.1,21.0,21.2,21.2,21.2,21.2,
1 21.2,21.2,21.2,21.3,21.3,21.3,21.3,21.1,21.2,21.2,21.3,21.3,21.1,
1 21.1,21.3,21.2,21.7,21.7,21.2,21.2,21.2,21.1,21.2,21.2,21.2,21.2,
1 21.5,21.3,21.5,21.2,21.3,21.3,21.2,21.2,21.3,21.1,21.2,21.2,21.2,
1 21.2,21.4,21.2,21.2,21.2,21.2,21.3,21.3,21.3,21.4,21.2,21.31.3/

```

N=0

```

C DO 200 I=1,NREC
C KREC=IREC+I-1
C PRINT *,IREC,KREC,N,NREC
C CALL DIRECT (KREC,CATFIL,2,I,IZO)
C DO 150 J=1,IFREQ
C N=N+1
C IZO (N)=ISPACE(J ,1)
150 CONTINUE

```

ORIGINAL PAGE IS
OF POOR QUALITY


```

300    CONTINUE
      NDZDT0=IDZDT0
      IF(IADMIN.LT.0) NDZDT0=IDZDT0-IADMIN
C      DO 100 I=1,2
C      SLOPE(I)=1.
C      INTRCP(I)=1.
100    CONTINUE
C      C=1.
C      BSA=1.
      DO 1000 I=1,NPEAKS
      ISTRT=OUTPUT(3,I)
      DO 400 J=ISTRT,2,-1
      IF(INPUT(J).GT.NDZDT0 ) GOTO 400
      TSTRT= J
      GOTO 500
400    CONTINUE
      TSTRT=0
C      HIT THE BEGINNING OF THE GROUP BUT STILL GOING DOWN; CHECK
C      IPTOLD FOR STARTING TIME
C      ISTRT=170
C425    CONTINUE
C      IF(ICOUNT.GT.1) THEN
C      DO 450 J=ISTRT,2,-1
C      IF(IPTOLD(J-1).LT.IPTOLD(J) ) GOTO 450
C      TSTRT=600.*(ICOUNT-1)-(170.-J+1)
C      GOTO 500
C450    CONTINUE
D      TYPE *, ' ERROR NO BOTTOM '
C      END IF
C      TSTRT=600.*(ICOUNT-1)
500    CONTINUE
      TIME=(OUTPUT(8,I )-TSTRT)
C      CALL CGL(37,3,0)
      CALL CGL(33,TSTRT,350.)
C      CALL CGL(37,5,0)
      CALL CGL(33,OUTPUT(8,I),350.)
      ZNOT=IZO( ISTRT ) *SLOPE(1)+INTRCP(1)
C      DO 550 J=1,NUM
C      IF(IDKREC.EQ.IEPC(J)) THEN
C      ZNOT=ZOEPC(J)
C      GOTO 575
C      ENDIF
550    CONTINUE
575    CONTINUE
      DZDT=OUTPUT(2,I)*SLOPE(2)+INTRCP(2)
      VNTVOL=C*TIME*(DZDT
                                     )/IFREQ/ZNOT**2
      DATCAR(1,I)=VNTVOL
      IF(I.EQ.1) THEN
      BPM=IFREQ*60./((OUTPUT(3,2)-OUTPUT(3,1)))
      ELSE
      BPM=IFREQ*60./((OUTPUT(3,I)-OUTPUT(3,I-1)))
      END IF
C      ELSE
C      END IF

```

```

CC=VNTVOL*BPM/1000.
CI=CO/BSA
C HI=0.0
  DATCAR(2,I)=CO
  DATCAR(3,I)=CI
  DATCAR(4,I)=DZDT
  DATCAR(5,I)=OUTPUT(3,I)

C
C   SAVE TIME FACTORS FOR HI FROM DZDT DATA AND DZDT PEAK VALUES
C   WHERE HI = DATCAR(4,I)/[Q-DATCAR(5,I)]
C   AND Q IS TIME OF Q WAVE BEFORE DZDT PEAK
C
D   WRITE(2,*      ) OUTPUT(2,I),TIME,VNTVOL,TSTRT,J,ISTRT,
D   1 ZNOT,C,DZDT,BPM,SLOPE,INTRCP,IZO(ISTRT)
D   WRITE(2,*      )(OUTPUT(3,I)-1.)/IFREQ,VNTVOL,CO,CI,HI,
D   1 TIME/IFREQ,DZDT,BPM,ZNOT,C
C   PRINT *,I,TIME,' ( ',OUTPUT(8,I), TSTRT,' ) ',ZNOT,
C   1 ' ( ',OUTPUT(3,I),' ) ',DZDT,' ( ',OUTPUT(2,I),' ) '
  600   CONTINUE
  1000  CONTINUE
C   OLDOUT=OUTPUT(3,NPEAKS)
C   DO 4000 J=431,600
C   IPTOLD(J-430)=INPUT(J)
C4000  CONTINUE
  RETURN
  END
  SUBROUTINE STRVOL (KPEAKS,DATSD,IPEAKS)
  PARAMETER (IFREQ=200)
  DIMENSION DATCAR(5,50),DATSD(4),IPEAKS(KPEAKS)
  COMMON /PLOT/X(1000),Y(1450)
  EQUIVALENCE (DATCAR,Y(1201))
  DO 1500 I=1,4
  SUM=0.0
  SUMSQ=0.0
  DO 1000 J=1,KPEAKS
  SUM=SUM+DATCAR(I,IPEAKS(J))
  SUMSQ=SUMSQ+DATCAR(I,IPEAKS(J))**2
1000  CONTINUE
  DATCAR(I,1)=SUM/KPEAKS
  DATSD(I)=SQRT(KPEAKS*SUMSQ-SUM*SUM)/KPEAKS
1500  CONTINUE
  RETURN
  END
  SUBROUTINE HICALC (NPEAKS,IFREQ,IECG,KREC,NREC,CATFIL)
  PARAMETER (JFREQ=200)
  VIRTUAL IECG(9000 )
  DIMENSION OUTPUT(10,70),DATCAR(5,50),INPUT(15*JFREQ)
  CHARACTER*14 CATFIL
  COMMON /PLOT/X(1000),Y(1450)
  EQUIVALENCE (Y(501),OUTPUT),(Y(1201),DATCAR),(X,INPUT)
  N=0
C   DO 200 I=1,NREC
C   IREC=KREC+I-1
C   TYPE *,CATFIL

```

ORIGINAL PAGE IS
OF POOR QUALITY

```

C      CALL DIRECT (IREC,DATFIL,4,I,IECG)
C      DO 150 J=1,IFREQ
C      N=N+1
C      IECC (N)=ISPACE( J      ,1)
150    CONTINUE
200    CONTINUE
      DO 1000 I=1,NPEAKS
      IEND=DATCAR(5,I)
      IF(I.NE.1) THEN
      ISTRT=(IEND+DATCAR(5,I-1))/2
      ELSE
      ISTRT=(1+IEND)/2
      END IF
      ICOUNT=IEND-ISTRT+1
      DO 500 J=1,ICOUNT
      INPUT(J)=IECG(ISTRT+J-1+6000)
500    CONTINUE
      ITIME=ISTRT+IMAX(INPUT,ICOUNT)-1
      TIME=ITIME
      CALL CCL(33,TIME,150.)
      TIME=(DATCAR(5,I)-TIME)/IFREQ
      IF(TIME) 600,600,700
600    CONTINUE
      DATCAR(4,I)=0.0
      GOTO 1000
700    CONTINUE
      DATCAR(4,I)=DATCAR(4,I)/TIME
1000   CONTINUE
      RETURN
      END
      FUNCTION IMAX(INPUT,ICOUNT)
      DIMENSION INPUT(ICOUNT)
      IMAX=ICOUNT
      DO 1000 I=ICOUNT-1,1,-1
      IF(INPUT(I).GT.INPUT(IMAX)) THEN
      IMAX=I
      END IF
1000   CONTINUE
      RETURN
      END

```



Report Documentation Page

1. Report No. NASA TM-101020		2. Government Accession No.		3. Recipient's Catalog No.	
4. Title and Subtitle A Computer Program for Processing Impedance Cardiographic Data: Improving Accuracy Through User-Interactive Software				5. Report Date October 1988	
				6. Performing Organization Code	
7. Author(s) Patricia S. Cowings, Karen Naifeh,* and Chet Thrasher*				8. Performing Organization Report No. A-88260	
				10. Work Unit No. 199-22-22-32	
9. Performing Organization Name and Address Ames Research Center Moffett Field, CA 94035 *University of California, San Francisco, CA 94122				11. Contract or Grant No.	
				13. Type of Report and Period Covered Technical Memorandum	
12. Sponsoring Agency Name and Address National Aeronautics and Space Administration Washington, DC 20546-0001				14. Sponsoring Agency Code	
15. Supplementary Notes Point of Contact: Patricia S. Cowings, Ames Research Center, MS 239A-2 Moffett Field, CA 94035 (415) 694-5724 or FTS 464-5724					
16. Abstract This report contains the source code and documentation for a computer program used to process impedance cardiography data. The cardiodynamic measures derived from impedance cardiography are ventricular stroke volume, cardiac output, cardiac index and Heather index. The program digitizes data collected from the Minnesota Impedance Cardiograph, electrocardiography (ECG), and respiratory cycles and then stores these data on hard disk. It computes the cardiodynamic functions using interactive graphics and stores the means and standard deviations of each 15-sec data epoch on floppy disk. This software was designed on a Digital PRO380 microcomputer and used version 2.0 of P/OS, with (minimally) a 4-channel 16 bit analog/digital (A/D) converter. Applications software is written in Fortran 77, and uses Digital's Pro-tool Kit Real Time Interface Library (PRTL), CORE Graphic Library (CGL), and laboratory routines. Source code can be readily modified to accommodate alternative detection, A/D Conversion and interactive graphics. The object code utilizing overlays and multitasking has a maximum of 50 Kbytes.					
17. Key Words (Suggested by Author(s)) Impedance cardiography, Cardiodynamics, Cardiac output, Cardiac index, Heather index, Stroke volume			18. Distribution Statement Unclassified-Unlimited Subject Category - 52		
19. Security Classif. (of this report) Unclassified		20. Security Classif. (of this page) Unclassified		21. No. of pages 57	
				22. Price A04	